



Company Headquarters (Baoding)

- Research, development design, laboratory
- Marketing, sales
- Quality control
- Manufacture
- Service
- Warehouse, logistics


Company History

- Founded in 1997
- Incorporate a joint venture company with Amerihua International Enterprises Inc. in 2001
- Pass the ISO 9001:2000 quality certification system in 2003
- Peristaltic pumps and syringe pumps obtain CE certifications in 2004
- Become a wholly foreign-owned enterprise in 2010

Company Mission

- Excellence because of profession: Longer only develops precision pumps
- Progress because of concept: Innovation and Precision
- Growth because of system: Perfect system results in success.

A Complete Peristaltic Pump




Peristaltic Pump Working Principle



Peristaltic Pump Features

**Non-contamination, easy to clean,
low-shear, low maintenance**

A Complete Syringe Pump



Syringe Pump Working Principle




Table of Contents


| | | |
|---|-----------------------------|---------|
| Peristaltic Pump Head | | 01 - 09 |
| Basic Peristaltic Pump | Flow Rates | |
| BQ50 - 1J | 0.0002mL/min - 20mL/min | 10 |
| BT100 - 2J | 0.0002mL/min - 380mL/min | 11 |
| BT300 - 2J | 0.07mL/min - 1140mL/min | 12 |
| BT600 - 2J | 0.07mL/min - 3000mL/min | 13 |
| WT600 - 2J | 4.2mL/min - 6000mL/min | 14 |
| WT600 - 3J | 4.2mL/min-6000mL/min | 15 |
| Flow Rates Peristaltic Pump | | |
| LEAD - 2 | 0.005mL/min - 380mL/min | 16 |
| BT100 - 1L | 0.002mL/min - 500mL/min | 17 |
| Dispensing Peristaltic Pump | | |
| BT100 - 1F | 0.0002mL/min - 500mL/min | 18 |
| BT300 - 1F | 0.07mL/min - 1140mL/min | 19 |
| WT600 - 1F | 0.7mL/min - 6000mL/min | 20 |
| Industrial Peristaltic Pump | | |
| YT600 - 1J | 600mL/min - 11000mL/min | 21 |
| WT600 - 4F | 100mL/min - 11000mL/min | 22 |
| Batch Transferring Pump | | |
| JL350 - 2J | ≤35 (L/min) | 23 |
| OEM | | 24 |
| Tubing | | 25-28 |
| Fitting | | 29 |
| Accessory | | 30 |
| Laboratorial Syringe Pump | | |
| TJ /TS Series Syringe Pump | 0.764nL/min - 1325 μ L/min | 31-34 |
| LSP01 - 1A(2A)、 LSP04 - 1A | 0.166nL/min - 54.15mL/min | 35 |
| LSP02(10) - 1B | 0.831nL/min - 150.5mL/min | 36 |
| LSP01 - 1C | 0.831nL/min - 10.84mL/min | 37 |
| LSP01 - 1BH | 0.09 μ L/min - 124.36mL/min | 38 |
| Micro Gear Pump | | |
| WT3000 - 1JA | 85.7mLmin - 2571.4mLmin | 39 |
| WT3000 - 1JB | 85.7mLmin - 2571.4mLmin | 40 |
| WT3000 - 1FA | 85.7mLmin - 2571.4mLmin | 41 |
| WT3000 - 1FB | 85.7mLmin - 2571.4mLmin | 42 |
| Diluters & Dispensers System | | |
| SMD02 - 1 | | 43 - 44 |
| SMD01 - 1 | | 43 - 44 |

NEW
Diluters & Dispensers
System



SMD02-1

Powerful



LSP01 - 1BH




NAVID KALA DARAN CO.

Easy-load Pump Head

YZ II 15(YZ II 25)


Patent No.:200620026529.1.

Change tubing easily and rapidly; Compact size. The tubing retention set is a linkage, It is convenient to load the tubing; Accept several tubing sizes for a wide flow range; The rollers adopt high quality materials and has good wearing property; The housing material is PESU. It has perfect rigidity and structure property. It is stable and autoclavable (200°C)



YZ1515x(YZ2515x)

Change tubing easily and rapidly; Automatic tubing retention; Suitable for several tubing sizes; The rollers adopt high quality materials and has good wearing property; The housing of the pump head has perfect rigidity and structure property, and its material is PESU which is stable and autoclavable (200 °C)



Tubing Setting



1 Press down the triggers according to the arrow heads directions and open the compression block



2 Load the tubing




3 Press the triggers according to the arrow heads directions and close the compression block


| Type | Tubing | Max. flow rates (mL/min) | Speed (rpm) | Housing Material | Part Number | | Weight (kg) |
|----------|-----------------------------------|-----------------------------|----------------|---------------------|--|------------------|-------------|
| | | | | | SS Rollers | RULON® J Rollers | |
| YZ1515x | 13" 14" 19" 16" 25" 17" 18" | 2200 | ≤ 600 | PESU | 0501521 (3 Rollers) 0501522 (6 Rollers) | | 0.4 |
| YZ2515x | 15" 24" | 1600 | | | 0501531 (3 Rollers) | | |
| YZ II 15 | 13" 14" 19" 16" 25" 17" 18" | 2200 | | | 0501552 | 0501551 | 0.35 |
| YZ II 25 | 15" 24" 35" 36" | 3000 | | | 0501562 | 0501561 | |



Micro Flow Rates Multi - channel Pump Head



DG-1, DG-2 Trigger Operation



Features



- Multi-channel transfer
- Change tubing easily
- Fix tubing easily
- Occlusion can be adjusted slightly by ratchet wheel to meet different tubing wall thickness requirements.
- 6-roller and 10-roller pump heads available.
- More rollers reduce pulsation and flow rates slightly.
- Skillful trigger design, convenient to open the cartridges (DG-1, DG-2)
- Rollers adopt 316 stainless steel which has chemical resistance to organic sdntion, acid and alkali.

Introduction

DG series pump heads are designed for small flow rates, multi-channel fluids transfer. Easy to change and fix the tubing. Occlusion can be adjusted slightly. The rollers adopt high quality materials. The pump head consists of base, rotor assembly and easily dismantled cartridge.

| Type | Tubing | Max. Flow Rates (mL/min) | | Speed (rpm) | Cartridge Material | Part Number | | Weight (kg) | |
|---------|---|--------------------------|---------------------|-------------|---------------------|-------------|------------|-------------|------------|
| | | 6 Rollers | 10 Rollers | | | 6 Rollers | 10 Rollers | 6 Rollers | 10 Rollers |
| DG - 1 | Inner Diameter ≤ 3.17 mm Wall Thickness 0.8 - 1 (mm) | 48 (per channel) | 32 (per channel) | ≤100 | POM 04011010001 | 0501101 | 0501102 | 0.21 | 0.20 |
| DG - 2 | | | | | | 0501111 | 0501112 | 0.27 | 0.26 |
| DG - 4 | | | | | | 0501121 | 0501122 | 0.40 | 0.39 |
| DG - 6 | | | | | | 0501131 | 0501132 | 0.54 | 0.51 |
| DG - 8 | | | | | PVDF 04011010002 | 0501141 | 0501142 | 0.67 | 0.63 |
| DG - 12 | | | | | | 0501151 | 0501152 | 0.95 | 0.88 |

Medium and Small Flow Rates Multi-channel Pump Head



Features

- Cartridges and base can be separated. It's convenient to load the tubing
- Suitable for several tubing sizes for a wide flow range
- Occlusion can be adjusted slightly to meet different tubing wall thickness requirements
- The material of the roller is stainless steel
- The material of the cartridge working surface is POM. It has perfect self lubricating property to reduce tubing wear abrasion

Pump Head Type



DG15 - 24



DG15 - 28



DG15 - 48

- DG15 - 24: 2 Channels 4 Rollers
- DG15 - 28: 2 Channels 8 Rollers
- DG15 - 48: 4 Channels 8 Rollers

Introduction

Patent No.: ZL 200520026479.2

DG15 series pump heads are designed for medium and small flow rates, multi-channel fluids transfer. Easy to change and fix the tubing. Occlusion can be adjusted slightly. The rollers adopt high quality materials. The pump head consists of base, rotor assembly and easy-to-load cartridges.

| Type | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Cartridge Material | Rollers Material | Rollers Quantity | Part Number | Weight (kg) |
|-----------|--|-----------------------------|----------------|-----------------------|---------------------|---------------------|-------------|-------------|
| DG15 - 24 | 16" 25" 17" | 1800 | ≤ 600 | POM | SS | 4 | 0501162 | 0.82 |
| DG15 - 28 | Inner Diameter ≤ 3.17mm Wall Thickness 0.8 - 1 (mm) | 75 | ≤ 100 | | | 8 | 0501172 | 0.67 |
| DG15 - 48 | | | | | | 13" 14" | 0501182 | 0.87 |


Standard Pump Head

Classical



Introduction

Stainless steel or high performance plastic rollers available
Transparent PC housing, visible operation



Tubing loading Procedure



A

Open the pump head housing. Press the tubing using the loading key



B

Rotate the loading key to load the tubing between rollers and housing



C

Close the pump head housing.

PERISTALTIC PUMP


Features

- Stable operation and stackable
- Compact size
- Dismounting the pump head from the drive is needed when loading the tubing
- Standard pump head, fixed occlusion
- Ideal for OEM

| Type | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Housing Material | Rollers Quantity | Part Number | | Weight (kg) |
|------|--------|--------------------------|-------------|------------------|------------------|-------------|------------------|-------------|
| | | | | | | SS Rollers | PET - TX Rollers | |
| BZ25 | 24" | 1600 | ≤ 600 | PC | 3 | 0501002 | 0501001 | 0.31 |

Low Pulse Pump Head



DMD25



Features

- Suitable for accurate dispensing
- Track assembly and base are independent. Easy to load tubing
- Accept many sizes of tubing to realize wide range flow rates transferring
- Twin offset track combing with Double-Y tube elements reduce the pulsation successfully. The accuracy can reach $\pm 0.5\%$
- Hay cutter style fixing prevents the tubing from quick abrasion caused by shock
- The material of track assembly is POM

Pulsation Reducing Sketch




Tubing changing procedure



Lift the two levers at both sides of the pump head and remove the track.

Loose the knobs, open the hay cutter style fixing by turning it through one quarter turn. Put the Double-Y tubing on the rollers and locate it over the mushroom pegs on both sides of the pump head.



Close the hay cutter style fixing. Tighten the knobs on both sides of the pump head.

Replace the track and secure it by closing the two levers.



Introduction

DMD25 pump head is designed for high accuracy dispensing. Twin offset track design and Double-Y tubing ensure the high dispensing accuracy. Hay cutter style fixing prevents the tubing from quick abrasion caused by shock. The pump head adopts stainless steel, aluminum and POM materials. The performance is reliable and stable. It is ideal for industry application.

Twin offset tracks reduce the pulsation, increase the dispensing accuracy.

| Type | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Track Material | Rollers Material | Rollers Quantity | Part Number | Weight (kg) |
|-------|--------|--------------------------|-------------|----------------|------------------|------------------|-------------|-------------|
| DMD25 | 119" | 150 | ≤ 350 | POM | SS | 6 | 0501051 | 2.5 |
| | 120" | 800 | | | | | 0501052 | |
| | 15" | 1500 | | | | | 0501053 | |
| | 24" | 2400 | | | | | | |
| | 35" | 2950 | | | | | | |
| | 36" | 4000 | | | | | | |

Quick-load Pump Head

KZ25

Flow Rates : 200 - 6000 mL/min



Adjusting Screw


Part Number

0501251 Long Shaft (for stack, connect drive directly)

0501252 Short Shaft (for stack)

Introduction

Stainless steel or RULON®J rollers available
PPS compression block. Self-lubrication to reduce tubing wear.
Good rigidity for stable size. Good chemical compatibility and high temperature resistance
High quality materials form a perfect combination. That offers excellent technical and service performance




Tubing loading Procedure



A

Rotate lever to loosen the compression block, and then lift it.



B

Load tubing on the rollers and the two V-shaped slots of tubing retainers. Put down the compression block.



C

Press the compression block and rotate lever to lock the compression block. Adjust the screws to fix the tubing.


Features

Cartridges and base can be separated. It's convenient to load the tubing

Pump heads can be stackable

Compact size and a wide flow range

Ideal for OEM



| Type | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Rollers Material | Housing Material | Compression Block Material | Rollers Quantity | Part Number | Weight (kg) |
|------|--------|--------------------------|-------------|------------------|------------------|----------------------------|------------------|-------------|-------------|
| KZ25 | 15" | 1800 | ≤ 600 | SS | PC | PPS | 3 | 0501251 | 0.79 |
| | 24" | 3500 | | RULON®J | | | | 0501252 | |
| | 35" | 5000 | | | | | | | |
| | 36" | 6000 | | | | | | | |

Industrial Easy-load Pump Head

YZ35-13



Tubing loading Procedure



A

Following the indication of arrow, rotate lever to left to open pump head.



B

Load the tubing



C


Following the indication of arrow, rotate lever to right to close pump head.



Knob

Rotate knob to micro-adjusted the gap between compression block and rotor .

Installation Drawing



Features

Pump heads can be stackable

Change tubing easily and rapidly

Automatic tubing retention

Suitable for several tubing sizes

The rollers adopt high quality materials and has good wearing property

The housing of the pump head has perfect rigidity and structure property

Its material is PSU which has stable performance and high temperature resistance (150 °C)

| Type | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Rollers Material | Rollers Quantity | Housing Material | Part Number | Weight (kg) |
|-----------|---------|--------------------------|-------------|------------------|------------------|------------------|--------------------|-------------|
| YZ35 - 13 | 73" 82" | 11000 | ≤ 600 | SS | 3 | PSU | 0501591 0501592 | 1.65 |

Industrial Quick-load Pump Head

KZ35



Tubing Fittings/Clamps Installation Procedure



A

Turn the left and right lever and take out the compression block.



B

Insert the tubing fittings/clamps with tubing into the slots on both sides.




C

Put the compression block on the tubing and turn the left and right levers to lock the compression block.

PUMP HEADS

Installation Drawing



Features

Key parts adopt 304 stainless steel. Suitable for pharmaceutical and food industry.

Compression block and base are separate. Easy to load the tubing.

Accept 73" and 82" tubing. The flow rate range is wide.

Two pump heads can be stacked.

The advantage of selecting tubing clamps: No residue fluids in the tubing

The advantage of selecting tubing fitting: Save tubing cost. The low cost tubing can be connected outside of the pump head.



| Pump Head | Tubing | Max. Flow Rates (mL/min) | Speed (rpm) | Tubing Clamp/Fitting Material | Pump Head Material | Roller Quantity | Part Number | Weight (kg) |
|-----------|---------|--------------------------|-------------|-------------------------------|---------------------|-----------------|--|-------------|
| KZ35 | 73" 82" | 11000 | <600 | PVDF | 304 Stainless Steel | 3 | 0501261 0501262 0501263 0501264 | 3.7 |


Other Pump Head



Part Number
0501451 (Black)
0501452 (White)

WX10-14 Patent No: ZL200420092094. 1


Features: Trigger, flip-open, small flow rates



Part Number
0501211

JY15-12


Features: Compact size, high flow rates, thicker wall thickness tubing are acceptable



Part Number
0501401 (Soft spring, with tubing connector)
0501402 (Hard spring, with out tubing connector)
0501403 (Soft spring, with tubing connector)
0501404 (Hard spring, with out tubing connector)

TH15

Features: Flexible rotor assembly, occlusion can be adjusted automatically, which prolong the tubing life



| Type | Tubing (mm) | Max. Flow Rates (mL/min) | Motor speed (rpm) | Rollers Quantity | Rollers Material | Part Number | Housing Material |
|-----------|--|--------------------------|-------------------|------------------|------------------|--|------------------|
| WX10 - 14 | Wall Thickness: 0.8 - 1.0 Inner Diameter: ≤ 3.17 | 24 | ≤ 60 | 4 | POM | 0501451 0501452 | ABS |
| JY15 - 12 | 25" 17" | 170 | ≤ 100 | 2 | PET | 0501211 | PPS |
| TH15 | 13" 14" 19" 16" 25" | 167 | ≤ 100 | 2 | POM PVDF | 0501401 0501402 0501403 0501404 | PPS |

Basic (Micro) Peristaltic Pump

Flow Rates




BQ50-1J



BQ50-1J peristaltic pump combines low flow and compact size into a single pump that mounts in the equipment rack or panel. The pump provides flow rates from 0.0002 to 20.0 mL/min. And it can be operated through hand-held remote controller, analog interface or RS485. White and black colors are available.

Features



OEM Installation Sketch

- Prime Function
 - Fast Filling
 - Fast Emptying
- Memory function
- Ideal for OEM
 - Compact size
 - Light weight
- Installation Modes
 - Inlaid
 - Rack
- Control Modes
 - Hand-held remote controller
 - Analog external control
 - Communication external control
- Tubing Loading
 - Flexible
 - Reliable



Specifications

| |
|--|
| Speed: 1 to 50 rpm, reversible |
| Speed precision: 1 rpm (External control 0.1 rpm) |
| Speed control: Hand-held remote controller |
| Display: Hand-held remote controller displays speed range and running status |
| External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, 4 to 20 mA and 0 to 10 kHz speed control |
| Communication interface: RS485 |
| Power supply: DC 12 V/1A adapter |
| Power consumption: < 10 W |
| Operating condition: Temperature 0 - 40 °C Relative humidity < 80% |
| Drive dimensions (L×W×H) : 135×72×72 (mm) |
| Controller dimensions (L×W×H) : 105×50×16 (mm) |
| Drive weight: 0.5 kg |
| IP rating: IP 31 |

| Standard Configuration | Optional Accessories |
|--------------------------------|--|
| 1. Pump Head | 1. Small V-base |
| 2. Drive | 2. Polished Stainless Steel Frame Tube |
| 3. Hand-held Remote Controller | 3. Fixed Plate |
| 4. Adapter | |
| 5. Data Wire | |




| Drive | Part Number | Pump Head | Flow Rates (mL/min) | Tubing (mm) | Weight (kg) |
|-----------|-------------|-----------|---------------------|--|-------------|
| BQ50 - 1J | 0502021 | WX10 - 14 | 0.0002 - 20 | Wall Thickness: 0.8 - 1.0 | 0.5 |
| | 0502022 | | | Inner Diameter: ≤ 3.17 ex. 0.5×0.8 1×1 2×1 2.4×0.8 3×1 etc. | |

Basic Peristaltic Pump

BT100-2J

Acceptable pump heads for BT100-2J are YZ1515x, YZ2515x, YZ1115, YZ1125, DG-1 and DG-2. The pump provides flow rates from 0.0002 to 380 mL/min. The speed can be adjusted manually or automatically through external control. Store the running parameters automatically. Easy to operate.



Applicable pump heads



YZ1515x
YZ2515x



YZ1115
YZ1125



DG-1



DG-2



Standard external control interface
—Order separately



Functions and Features

Applicable pump heads: YZ1515x, YZ2515x, YZ1115, YZ1125, DG - 1, DG - 2

Max key: For fast filling and emptying

Communication function: RS485

External control: Start/stop and cw/ccw control, 0 - 5 V, 0 - 10 V, 4 - 20 mA, 0 - 10 kHz speed control

Adjust speed manually or automatically through external control interface

Memory function: Storing the running parameter automatically

Membrane keypad, easy to operate

Good performance and economical price

At normal lab conditions

Specifications


- Speed: 0.1 to 100 rpm, reversible
- Speed precision: 0.1 rpm
- Speed control: Membrane keypad
- Display: 3-digit LED displays current rpm
- External control: Start/stop, cw/ccw control, and 0 - 5 V, 0 - 10 V, 4 - 20 mA and 0 - 10kHz speed control
- Power supply: AC 90 V - 260 V 50/60 Hz
- Power consumption: ≤ 30 W
- Operating condition: Temperature 0 to 40 °C
Relative humidity < 80%
- Dimensions (L×W×H) : 232×142×149 (mm)
- Drive weight: 2.3 kg
- IP rating: IP 31

| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|---------------------|--|-------------------------------|-------------|
| BT100 - 2J | 0502111 | YZ1515x、YZ 1115 | 13# 14# 19# 16# 25# 17#18# | 0.07 - 380 | 3.0 |
| | | YZ 2515x | 15# | 1.7 - 150 | 3.0 |
| | | YZ 1125 | 15# | 1.7 - 150 | 3.0 |
| | | DG - 1 (6 Rollers) | | 0.00025 - 48 | 2.5 |
| | | DG - 1 (10 Rollers) | Wall Thickness: 0.8 - 1.0 (mm) Inner Diameter: ≤ 3.17 mm | 0.0002 - 32 | 2.5 |
| | | DG - 2 (6 Rollers) | | 0.00025 - 48 (Single channel) | 2.6 |
| | | DG - 2 (10 Rollers) | | 0.0002 - 32 (Single channel) | 2.6 |

Basic Peristaltic Pump

BT300-2J

Acceptable pump heads for BT300-2J are YZ1515x, YZ2515x, YZII15, YZII25. The pump delivers flow rates from 0.07 to 1140 mL/min. The speed can be adjusted manually or automatically through external control interface. Store the running parameters automatically. Easy to operate.



Applicable Pump Heads



YZ1515x
YZ2515x



YZII15
YZII25

PERISTALTIC PUMP

Specifications

| |
|--|
| Speed: 1 to 300 rpm, reversible |
| Speed precision: 1 rpm |
| Speed control: Rotary encoded switch |
| Display: 3-digit LED displays current speed |
| External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, 4 to 20 mA and 0-10 kHz speed control |
| Power supply: AC 90 V - 260 V 50/60 Hz |
| Power consumption: < 48 W |
| Operating condition: Temperature 0 to 40 °C Relative humidity < 80% |
| Dimensions (L W H): 285 207 180 (mm) |
| Drive weight: 3.6 kg |
| IP rating: IP 31 |

Functions and Features


- Applicable pump heads: YZ1515x, YZ2515x, YZII15, YZII25
- Prime button for fast filling and emptying
- Adjust speed manually or automatically through external control
- Memory function, storing the running parameters automatically
- Equipped with Longer dispensing controller for filling function
- At normal lab conditions



| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|------------------|-----------------------------|---------------------|-------------|
| BT300 - 2J | 0502311 | YZ1515x, YZII 15 | 13# 14# 19# 16# 25# 17# 18# | 0.07 - 1140 | 4.0 |
| | | YZ2515x, YZII 25 | 15# 24 # | 1.7 - 870 | 4.0 |

Basic Peristaltic Pump

Flow Rates



BT600-2J

Applicable pump heads for BT600-2J are YZ1515x, YZ2515x, YZ1115, YZ1125. The pump delivers flow rates from 0.07 to 3000 mL/min. The speed can be adjusted manually or automatically through external control. Easy to operate. Dispensing function can be realized by connecting Longer dispensing controller FK-1C.



Applicable Pump Heads



YZ1515x
YZ2515x




YZ1115
YZ1125

Specifications

| |
|---|
| Speed: 1 to 600 rpm, reversible |
| Speed precision: 1 rpm |
| Speed control: Rotary encoded switch |
| Display: 3 - digit LED displays current rpm |
| External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, 4 to 20 mA and 0-10 kHz speed control |
| Communication interface: RS485 |
| Power: AC 90 - 260 V 50/60 Hz |
| Power consumption: < 100 W |
| Operating condition: Temperature 0 to 40°C |
| Relative humidity < 80% |
| Dimensions (L×W×H) : 285×207×180 (mm) |
| Drive Weight: 3.8 kg |
| IP rating: IP 31 |

Functions and Features

- Acceptable pump heads: YZ1515x, YZ2515x, YZ1115, YZ1125
- Prime button: For fast filling and emptying
- Memory function, storing the running parameters automatically
- The speed can be controlled by external control analog input or pulse input
- Good performance and economical price
- Indoor use only
- Realize computer control through RS-485 communication interface.



| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|-----------------|----------------------------|---------------------|-------------|
| BT600 - 2J | 0502621 | YZ1515x, YZ1115 | 13# 14# 19# 16# 25# 17#18# | 0.07 - 2200 | 4.2 |
| | | YZ2515x | 15# 24# | 1.6 - 1600 | 4.2 |
| | | YZ1125 | 15# 24# 35# 36# | 1.7 - 3000 | 4.2 |

Basic Peristaltic Pump


WT600 - 2J

WT600-2J peristaltic pump adopts brushless DC motor. It has features of high output torque, low vibration, high efficiency and free maintenance. It can drive double pump heads. It accepts many pump heads and delivers flow rates from 4.2 to 6000 mL/min. The speed can be adjusted manually or automatically through external control.



Functions and Features

- Prime button: For fast filling and emptying
- Memory function: Store the speed and the address of the pump
- Communication function: Realize the computer control through RS485 interface
- Display: 3 digits LED displays the speed of the pump
- External control input: Control the speed, start/stop and direction of the pump
- Suitable for different pump heads and tubing for multiple application
- High output torque which can drive double pump heads
- Brushless DC motor, high output torque and free maintenance
- Operation mode: Switch, button and knob. Easy to operate



Applicable Pump Heads



YZ1515x
YZ2515x



DMD25



YZII15
YZII25



KZ25



DG15-24

Specifications

| |
|---|
| Speed: 60 to 600 rpm |
| Speed accuracy: $\pm 1.0\%$ |
| Speed resolution: 1.0 rpm |
| Output torque: $\geq 1.60 \text{ N} \cdot \text{M}$ |
| External control input: Start/stop and cw/ccw control, 0.5 to 5 V, 1 to 10 V, 4 to 20 mA and 1-10 kHz speed control |
| Communication interface: RS485 |
| Power supply: AC 176 - 264V 50Hz/60Hz |
| Power consumption: $\leq 200 \text{ W}$ |
| Operating condition: Temperature 0 to 40 °C |
| Relative humidity < 80% |
| Dimensions (L×W×H) : 285×207×180 (mm) |
| Drive weight: 5.2 kg |
| IP rating: IP 31 |

| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|----------------------------------|-----------------------------|---------------------------|-------------|
| WT600 - 2J | 0502631 | (1, 2) × YZ1515x | 13# 14# 19# 16# 25# 17# 18# | Single channel 4.2 - 2200 | 5.6 - 6.0 |
| | | (1, 2) × YZII 15 | | | |
| | | (1, 2) × YZ2515x | 15# 24# | Single channel 100 - 1600 | |
| | | (1, 2) × YZII 25 | 15# 24# 35# 36# | Single channel 100 - 3000 | 6.0 |
| | | KZ25 | | | |
| | | DG15 - 24 | 16# 25# 17# | 50 - 1800 | 6.0 |
| | | DMD25 ($\leq 350 \text{ rpm}$) | 119# 120# 15# 24# 35# 36# | 30 - 4000 | 7.7 |

* Max. speed for DMD25 is 350 rpm

Basic Peristaltic Pump

WT600 - 3J


WT600-3J peristaltic pump is high IP rating and high efficiency. It adopts brushless DC motor and speed closed loop control. It has features of high output torque, low vibration and free maintenance. It accepts multiple two channels pump heads and delivers flow rates from 4.2 to 6000 mL/min (single channel) .



Part Number
0502641 (0.5-5V)
0502642 (4-20mA)
0502643 (1-10V)
0502644 (1-10kHz)

Functions and Features


- High IP rating: Suitable for moist and dust working condition
- MAX key: For fast filling and emptying
- Memory function: Storing the running parameters automatically
- Communication function: Realize the computer control through RS485 interface
- Display: 3 digits LED displays the speed of the pump
- External control input: Control the speed, start/stop and direction of the pump
- Suitable for different pump heads and tubing
- High output torque which can drive double pump heads
- Brushless DC motor, high output torque, free maintenance
- Operation mode: Switch and membrane keypad. Easy to operate



Applicable pump heads



YZ1515x
YZ2515x



DMD25



YZII15
YZII25



KZ25



DG15-24


Specification

| |
|--|
| Speed: 60 - 600 (rpm) |
| Speed accuracy: $\pm 1.0\%$ |
| Speed resolution: 1.0 rpm |
| Output torque: $\geq 1.50 \text{ N}\cdot\text{M}$ |
| Display: 3 digits LED displays the speed of the pump; 4 pcs LED indicators display the working state |
| Memory function: Storing the running parameters automatically when the parameters don't change after 3 seconds |
| External speed control: 0.5 - 5 V, 1 - 10 V, 4 - 20 mA or 1 - 10 kHz are corresponding to 60 - 600 rpm |
| Communication interface: RS485 |
| Power supply: AC 176 - 264V 50Hz/60Hz |
| Power consumption: $\leq 190 \text{ W}$ |
| Operating condition: Temperature 0 to 40°C Relative humidity $< 80\%$ |
| Drive dimensions (L×W×H) : 273×190×272 (mm) |
| Drive weights: 7.3 kg |
| IP rating: IP55 |

| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|--|----------------------------------|-----------------------------|---------------------------|-------------|
| WT600 - 3J | 0502641 0502642 0502643 0502644 | (1,2) YZ1515x | 13# 14# 19# 16# 25# 17# 18# | Single channel 4.2 - 2200 | 7.7 - 8.1 |
| | | (1,2) YZII15 | | | |
| | | (1,2) YZ2515x | 15# 24# | Single channel 100 - 1600 | |
| | | (1,2) YZII25 | 15# 24# 35# 36# | Single channel 100 - 3000 | 8.1 |
| | | KZ25 | | 200 - 6000 | |
| | | DG15 - 24 | 16# 25# 17# | 50 - 1800 | |
| | | DMD25 ($\leq 350 \text{ rpm}$) | 119# 120# 15# 24# 35# 36# | 30 - 4000 | 9.8 |

* Max. speed for DMD25 is 350 rpm

Flow Rates Peristaltic Pump




LEAD-2

The applicable pump heads for Lead-2 are Lead15-48, Lead15-88, Lead15-24, Lead15-44 and max. 8 channels are available. The pump delivers flow rates from 0.005 to 380 mL/min. It has flow rates display and flow rates calibration functions. It also has functions of speed adjustment, start/stop, cw/ccw, prime (rinsing), tubing selection, RS485 communication, analog control, etc.



Applicable Pump Heads



LEAD15-24
(2 channels, 4 rollers)



LEAD15-44
(4 channels, 4 rollers)



LEAD15-48
(4 channels, 8 rollers)



LEAD15-88
(8 channels, 8 rollers)

Specifications

Speed: 1.0 to 100.0 rpm, reversible

Speed precision: 0.1 rpm

Speed control: Rotary encoded switch coupled with membrane keypad

Display: 128×32 graphic LCD displays current running status

External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, 4 to 20 mA and 0-10 kHz speed control

Communication: RS485

Power supply: AC 90 - 260V 50/60 Hz

Power consumption: < 50W

Operating condition: Temperature 0 to 40 °C
Relative humidity < 80%

Dimensions (L × W × H): 190×162×275 (mm)

Drive weight: 3.5 kg

IP rating: IP 31

Functions and Features

Vertical structure

Accept 4 special pump heads

Display can switch between flow rates and speed

Calibration function for more accurate flow rates

Computer control available through RS485 communication

Prime key for fast filling and emptying



| Drive | Part Number | Pump Head | Tubing | Flow Rates (Single channel:mL/min) | Weight (kg) |
|----------|-------------|-------------|---|---------------------------------------|-------------|
| LEAD - 2 | 0502197 | LEAD15 - 24 | 16# 25# 17# | 1 - 380 | 4.36 |
| | | LEAD15 - 44 | | 1 - 380 | 4.42 |
| | | LEAD15 - 48 | Inner Diameter ≤ 3.17 mm Wall Thickness 0.8 - 1.0 mm 13# 14# | 0.005 - 75 0.07 - 20 | 4.48 |
| | | LEAD15 - 88 | Inner Diameter ≤ 3.17 mm Wall Thickness 0.8 - 1.0 mm 13# 14# | 0.005 - 75 0.07 - 20 | 4.48 |

Flow Rates Peristaltic Pump


BT100-1L

The max. output torque of BT100-1L drive is 1.65 N·M. Multiple pump heads are available. The pump delivers flow rates from 0.002 to 500 mL/min. 128 × 64 graphic LCD displays running menus and parameter setting menus. The flow rates and the speed (rpm) can be displayed in the same screen. The display interfaces are friendly. It has flow rates calibration function to make the flow rates more accurate. It has many control modes. The external control can be realized through standard external interface or RS485 communication interface.



Functions and Features

Display can switch between flow rates and rpm
Calibration capability for more accurate flow rates
Graphic LCD together with rotary encoded switch make the pump easy to use
Higher torque to accept more pump heads and more channels
RS485 communication interface available. Offering communication protocol to customer and customizing software according to customer's requirements



Flow Rates Display High Torque

Applicable Pump Heads



DG-2
I
DG-24



DG15-24
DG15-28
DG15-48



YZ1515x
YZ2515x



YZII15
YZII25

Specifications

Speed: 1.0 to 100.0 rpm, reversible
Speed precision: 0.1 rpm
Speed control: Rotary encoded switch coupled with membrane keypad
Display: 128×64 graphic LCD displays current running status
External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, and 4 to 20 mA speed control
Communication interface: RS485
Power supply:
AC 220 V ± 10% 50Hz/60 Hz (standard)
AC 110 V ± 10% 50Hz/60 Hz (optional)
Power consumption: < 50 W
Operating condition: Temperature 0 to 40 °C
Relative humidity < 80%
Dimensions (L × W × H): 202×160×239 (mm)
Drive weight: 5.34 kg
IP rating: IP 31



| Drive | Part Number | Pump Head | Tubing | Flow Rates (Single channel: mL/min) | Weight (kg) |
|------------|-------------|---|---|--|-------------|
| BT100 - 1L | 0502122 | DG15 - 24 | 16# 25# 17# | 1 - 380 | 6.16 |
| | | DG15 - 28, DG15 - 48 | Wall Thickness: 0.8 - 1.0 (mm) Inner Diameter: ≤ 3.17 mm 13# 14# | 0.005 - 75 | 6.01 |
| | | DG - (2,4,6,8,12,16,24) (6,10 rollers) | Wall Thickness: 0.8 - 1.0 (mm) Inner Diameter: ≤ 3.17 mm | 0.0025 - 48 (6 rollers) 0.002 - 32 (10 rollers) | 5.82 - 8.70 |
| | 0502125 | (1,2,3,4) YZ1515x (1,2,3,4) YZ II 15 | 13# 14# 19# 16# 25# 17# 18# | 0.07 - 380 | 5.76 - 7.02 |
| | | (1,2,3,4) YZ2515x | 15# 24# | 1.7 - 270 | 5.76 - 7.02 |
| | | (1,2) YZ II 25 | 15# 24# 35# 36# | 1.7 - 500 | 5.70 - 6.04 |

Dispensing Peristaltic Pump


BT100-1F

Acceptable pump heads for BT100-1F are YZ1515x, YZ2515x, YZII15, YZII25 and DG-1, DG-2, DG-4. The pump delivers flow rates from 0.2 μ L to 500 mL/min and dispensing volume is from 0.01 mL to 9900 mL. Combining 128 \times 32 LCD display with membrane keypad and rotary encoded switch makes the operation easy and prompt. The speed can be adjusted manually or automatically through external control interface.



Functions and Features

| |
|--|
| Acceptable pump heads: YZ1515x, YZ2515x, YZII15, YZII25 and DG - 1, DG - 2, DG - 4 |
| Micro flow rates transferring: The smallest flow rate can reach 0.2 μ L/min |
| Back suction function: In dispensing mode, the pump runs reversely after finishing one filling operation to prevent the liquid from dripping |
| Operating mode: Membrane keypad and rotary encoded switch |
| Speed control: The speed can be adjusted manually or automatically through external control interface |
| Display: LCD displays current running status |
| Prime key: For fast filling and emptying |
| Memory function: Storing the running parameters automatically |
| Calibration function: The flow rates and the dispensing volume can be calibrated to increase the accuracy |
| Communication function: Realize computer control through RS485 communication interface |



Applicable Pump Heads



YZ1515x
YZ2515x



DG-4



YZII15
YZII25



DG-1
DG-2

Specifications

| |
|---|
| Flow rates: 0.2 μ L / min to 500 mL / min |
| Dispensing volume: 0.01 mL to 9.99 L |
| Speed: 0.1 rpm to 100.0 rpm |
| Copy number: 0 to 9999 |
| Pause time: 0.1 s to 99.9 min |
| Back suction time: 0 to 99.9 s |
| External control: start/stop control, direction control, flow rates control (0 - 5V, 0 - 10V, 4 - 20 mA, 0 - 10 kHz optional) |
| Signal output: start/stop, direction output and 0 - 10 kHz frequency output (corresponding to 0 to 100.0 rpm) with insulating function |
| Communication: RS485 |
| Power supply: AC 100 - 240V \pm 10% 50Hz/60Hz |
| Power consumption: < 40 W |
| Operating condition: Temperature 0 to 40°C Relative humidity < 80% |
| Dimensions (L \times W \times H): 285 \times 207 \times 180 (mm) |
| Drive weight: 3.8 kg |
| IP rating: IP 31 |



| Drive | Part Number | Pump Head | Tubing | Flow Rates (*/min) | Weight (kg) |
|------------|-------------|--|---|----------------------------------|-------------|
| BT100 - 1F | 0502132 | YZ1515x, YZ II 15 | 13" 14" 19" 16" 25" 17" 18" | 7 μL - 380 mL | 4.2 |
| | | YZ2515x | 15" 24" | 0.17 mL - 500 mL | |
| | | YZ II 25 | 15" 24" 35" 36" | | |
| | | DG - 1, DG - 2, DG - 4 (6 rollers) | Wall Thickness: 0.8 - 1.0 (mm) Inner Diameter: ≤ 3.17 mm | 0.25 μL - 48 mL (single channel) | 4.0 - 4.2 |
| | | DG - 1, DG - 2, DG - 4 (10 rollers) | | 0.2 μL - 32 mL (single channel) | |

Dispensing Peristaltic Pump

BT300-1F

Applicable pump heads for BT300-1F are YZ1515x, YZ2515x, YZ1115, YZ1125. The pump delivers flow rates from 0.07 to 1140 mL/min. 128 x 64 graphic LCD displays running menus and parameters setting menus. The display interfaces are friendly. It has many control modes. The external control can be realized through standard external control interface or RS485 communication interface.



Functions and Features

Acceptable pump heads: YZ1515x, YZ2515x, YZ1115, YZ1125

Set pause time to realize dispensing and filling automatically

Display can switch between flow rates and rpm


Calibration function to acquire more accurate flow rates

RS485 communication interface available. Offering communication protocol to customer and customizing software according to customer's requirements.



| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|----------------|-----------------------------|---------------------|-------------|
| BT300 - 1F | 0502322 | YZ1515x YZ1115 | 13# 14# 19# 16# 25# 17# 18# | 0.07 - 1140 | 4.7 |
| | 0502325 | YZ2515x YZ1125 | 15# 24# | 1.7 - 870 | 4.7 |

Flow Rates



Applicable Pump Heads



YZ1515x
YZ2515x



YZ1115
YZ1125

Specifications

Speed: 1 to 300 rpm, reversible

Speed precision: 1 rpm

Speed control: Rotary encoded switch coupled with membrane keypad

Display: 128 × 64 graphic LCD displays flow rates, speed or dispensing mode

Copy number: 1 to 999

Dispense volume: 0.1 mL to 99.9 L

Pause time: 1 to 999 s

Back suction angle: 0 to 360 18° increments

External control: Start/stop and cw/ccw control, 0 to 5 V, 0 to 10 V, and 4 to 20 mA speed control

Communication interface: RS485

Power supply:

AC 220 V ± 10% 50/60 Hz (standard)

AC 110 V ± 10% 50/60 Hz (optional)

Power consumption: < 50 W

Operating condition: Temperature 0 to 40 °C
Relative humidity < 80%

Dimensions (L × W × H): 202 × 160 × 239 (mm)

Drive weight: 4.3 kg

IP rating: IP 31

Dispensing Peristaltic Pump

WT600-1F

Acceptable various pump heads for WT600-1F are YZ1515x, YZ2515x, YZ1115, YZ1125, KZ25, DMD25 and DG15-24. The pump delivers flow rates from 0.7 mL to 6000 mL/min and dispensing volume is from 0.1 mL to 99.9 L. Combining 128 × 32 LCD display with membrane keypad and rotary encoded switch makes the operation easy and prompt. The speed can be adjusted manually or automatically through external control interface.



Functions and Features

Acceptable pump heads: YZ1515x, YZ2515x, YZ1115, YZ1125, DG15 - 24, KZ25, DMD25

Back suction function: In dispensing mode, the pump runs reversely after finishing one filling operation to prevent the liquid from dripping

Operating mode: Membrane keypad and rotary encoded switch

Speed control: The speed can be adjusted manually or automatically through external control

Display: LCD displays current running status and parameters

Prime key: For fast filling and emptying

Memory function: Storing the running parameters automatically

Calibration function: The flow rates and the dispensing volume can be calibrated to increase the accuracy


Flow Rates



Applicable Pump Heads



YZ1515x
YZ2515x



DMD25



YZ1115
YZ1125



KZ25



DG15-24

Specifications

Flow rates: 0.7 to 6000 mL / min

Dispensing volume: 0.1 mL to 99.9 L

Speed: 10 to 600 rpm

Copy number: 0 to 9999

Pause time : 0.1 s to 99.9 min

Back suction: 0 to 9.9 revolution

Flow rates calibrating time: 0.5 - 30 min

External control: start/stop control, direction control, flow rates control (0 - 5V, 0 - 10V, 4 - 20 mA, 0 - 10kHz optional)

Output interface: start/stop, direction output and 0 - 10 kHz frequency output

Communication: RS485

Power supply: AC 176V - 264V, 50Hz/60Hz

Power consumption: < 140 W

Operating condition: Temperature 0 to 40°C
Relative humidity < 80%

Drive Dimensions(L×W×H):

285×207×180 (mm)

Drive weight: 5.2 kg

IP rating: IP 31



| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|-----------------|-----------------------------|---------------------------|-------------|
| WT600 - 1F | 0502661 | (1,2) ×YZ1515x | 13# 14# 19# 16# 25# 17# 18# | Single channel 0.7 - 2200 | 5.6 - 6.0 |
| | | (1,2) ×YZ1115 | | | |
| | | (1,2) ×YZ2515x | 15# 24# | Single channel 17 - 1600 | |
| | | (1,2) ×YZ1125 | 15# 24# 35# 36# | Single channel 17 - 3000 | 6.0 |
| | | KZ25 | | | |
| | | DG15 - 24 | 16# 25# 17# | 8.2 - 1800 | 6.0 |
| | | DMD25 (≤350rpm) | 119# 120# 15# 24# 35# 36# | 5 - 4000 | 7.7 |

* Max. speed for DMD25 is 350 rpm


Tel: 86 - 312 - 3110087

Fax: 86 - 312 - 3168553

Internet
www.longerpump.com

E - mail: longer@longerpump.com

Industrial Peristaltic Pump



YT600-1J



Applicable pump heads



YZ35-13



KZ35



The pump delivers flow rates from 600 to 11000 mL/min. The speed can be adjusted manually or automatically through external control interface. Driven by DC motor YT600-1J has higher torque and can drive double pump heads. Suitable for industrial applications which need high pressure and flow rates.

Specifications

| | |
|---------------------|---|
| Speed | 60–600 (rpm) , reversible |
| Speed control | 10 turn potentiometer |
| Power supply | AC 220 V±10% 50/60 Hz |
| Power consumption | <400 W |
| Operating condition | Temperature 0 to 40°C Relative humidity <80% |
| External control | Start/stop control, speed control (1 - 10 kHz, 1 - 10 V, 4 - 20 mA optional, 4 - 20 mA is standard configuration) |
| Dimensions (L×W×H) | 325×236×193 (mm) |
| IP rating | IP 54 |


Functions and Features

- Acceptable pump heads: YZ35-13, KZ35
- Suitable for industrial applications
- Can be equipped with Longer Dispensing Controller for filling function when connecting single channel pump head

| Drive | Part Number | Pump Head | Tubing | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|-------------|---------|------------------------------|-------------|
| YT600 - 1J | 0502521 | YZ35 - 13 | 73" 82" | 600 - 11000 | 20 |
| | 0502522 | KZ35 | | 600 - 11000 | |
| | 0502523 | 2×YZ35 - 13 | | (Single Channel) 600 - 11000 | |
| | | 2×KZ35 | | (Single Channel) 600 - 11000 | |

Industrial Peristaltic Pump

WT600-4F



WT600-4F is a high IP rating and high efficient pump with high-power DC brushless motor, it is suitable for working in dust and damp industrial environment. It delivers flow rates from 100 to 11000 ml/min and dispensing volume from 100ml to 999 liters when loading YZ35 or KZ35 pump head. WT600-4F adopts 128 32 LCD to display all the running information and parameters and membrane keypad to set parameters. It is easy to operate. This pump can be controlled adopts manually or externally through standard external control module and RS485.



Specifications

Flow rates: 100 to 11000 ml/min
Dispensing volume: 100 ml to 999 liters
Speed: 10 to 600 rpm
Copy Number: 0 to 9999, "0" means running continuously
Pause time: 1s to 99.9 min
Back suction: 0 to 9.9 revolutions, (increment is 0.1 revolution)
Calibration time: 30 to 1800 seconds
External control input: Start/stop control, direction control, flow rates control (0 - 5V, 0 - 10V, 4 - 20 mA, 0 - 10kHz optional)
External control output: Start/stop output, direction output, 0 - 7.5 kHz frequency output corresponding to 0 to 600 rpm.
Communication: RS485
Power: AC 176 - 264 V, 50/60 Hz
Power consumption: < 300 W
Operating condition: Temperature 0 to 40 C, Relative humidity < 100%
Dimensions (L × W × H): 292×185×180 (mm)
Drive weight: 12 kg
IP Rating: IP65






Functions and Features

Acceptable pump head: YZ35-13, KZ35.
 Back suction: Prevent dripping of liquid when dispensing a certain liquid volume.
 Operating mode: Membrane keypad to set parameters and operation.
 Control mode: Pump can be controlled internally through membrane keypad or externally through external control modules.
 Display: 128 32 LCD displays all the running information and parameters.
 Prime function: Fast filling or emptying the tubing.
 Memory function: Save the parameters and working status automatically.
 Calibration: The flow rates and dispensing volume can be calibrated to improve the accuracy.

| Drive | Part Number | Pump Head | Housing Material | Tubing | Flow Rates (mL/min) | Weight (kg) |
|----------|-------------|-----------------|---------------------|---------|---------------------------------|-------------|
| WT600-4F | 0502662 | YZ35 - 13 (1,2) | PSU | 73# 82# | 100 - 11000 (Single Channel) | 13.65 |
| | 0502664 | KZ35 (1,2) | 304 Stainless Steel | | | 15.7 |
| | 0502663 | | | | | |
| | 0502661 | | | | | |

Batch Transferring Pump

JL350-2J



INDUSTRIAL

Specifications

Speed:

30 to 350 rpm, reversible

Speed precision: 0.6 rpm

Power supply:

AC 220 V \pm 10% 50/60 Hz


Power consumption: < 400 W

Operating condition: Temperature 0 to 40 °C
Relative humidity < 80%


Dimensions (L×W×H): 417×401×321 (mm)

Drive weight: 32 kg

IP rating: IP 55



Flow Rates



Tubing Loading



- A** 1. Turn the knob ① ② counter clockwise.
2. Turn the hand wheel to loose the compression block
3. Pull the knob ③ and take out the front cover



- B** 4. Select suitable tubing slots. Place the tubing between the roller and the compression block, then put the tubing in the corresponding tubing slots



- C** 5. Pull the knob ④, insert the front cover, tighten the knob ①
6. Turn the hand wheel and adjust the occlusion between the roller and the compression block
7. Tighten the knob ② and fix the compression block

Functions and Features

- High-power pumps, and the maximum flow rates can reach 35L/min
- The structure of concave working surface and convex roller realizes the tubing self-position function. It reduces abrasion and prolongs tubing life.
- Variable frequency speed control to AC motor with creeper gear
- Adjustable occlusion
- Strong drive
- Pump head operation is visible through transparent front cover

Tubing Fitting Structure



| Drive | Part Number | Pump Head | Tubing | Flow Rates (L/min) | Weight (kg) |
|------------|-------------|-----------|--------|--------------------|-------------|
| JL350 - 2J | 0502352 | KZ48/63 | 86" | 0.8 - 7.5 | 32 |
| | | | 88" | 1.3 - 15 | |
| | | | 90" | 2.7 - 32 | |
| | | | 92" | 3 - 35 | |

OEM Products

Longer company has accumulated rich experience in research and manufacture of peristaltic pump. Utilizing present resources which include software, hardware, materials and market, etc. we can provide reasonable and economical OEM products with professional design. Longer company also can design or manufacture special products to meet customers' requirements.






Utilize present products — Stable functions and short delivery period (preferred)

Custom design — To meet your special requirements (Batch purchasing)


OEM Applications

- Fluids sampling of inspection or testing instruments
- Feeding of fluids equipments
- Dispensing and filling fluids equipments
- Other fluids transfer

OEM Examples

| Type | Introduction | Functions and Features | Specifications |
|--|--|--|--|
|  OEM10/TH15 | OEM10 has four kinds of optional fixed speed which can deliver max. flow rates of 167 ml/min, it is one of ideal OEM product and has been widely used in environment monitoring, ferment, and filling industry, etc. | OEM10 adopts AC synchronous motor which has stable speed and high repeatability. Pump head with two spring rollers can reduce the abrasion of tubing and prolong the tubing life. | Speed: Four kinds of fixed speed of 5, 20, 60, 110 rpm, reversible Power supply: AC220V Power consumption: 15W Dimensions: 105.5 110 80 (mm) Weight: 0.7kg |
|  OEMBJ60-01/WX10 | The pump delivers flow rates from 0.001 to 24 mL/min. The speed can be adjusted by BCD DIP coding switch which is divided into 15 levels, or controlled by external pulse. | Suitable for tubing which the inner diameter is less than 3.17 mm and wall thickness is from 0.8 to 1.0 mm. Reasonable and fine shape, compact size, DC power supply, ideal for OEM. | Speed: 0 - 60 rpm, reversible Speed precision: 1 rpm Speed control: BCD switch, the speed is 1, 3, 5, 7 rpm when the speed is less than 10 rpm and the increment is 5 rpm when the speed is more than 10 rpm Power supply: 12V/1A DC Power consumption: Less than 10 W Dimension: 116 60 73 (mm) Weight: 0.55 kg |
|  OEMBJ100-01/JY15-12 | The pump delivers flow rates from 0.7 to 170 mL/min. The speed can be adjusted by BCD DIP coding switch which is divided into 15 levels, or controlled by external pulse. | Suitable for tubing with wall thickness of 1.6 mm. Low speed, high flow rates, long life. Reasonable and fine shape, compact size, DC power supply, ideal for OEM. | Speed: 0 - 60 rpm, reversible Speed precision: 1 rpm Speed control: BCD switch, the increment is 5 rpm when the speed is less than 50 rpm and the increment is 10 rpm when the speed is more than 50 rpm Power supply: 12V/2A DC Power consumption: Less than 25 W Dimension: 118 60 86 (mm) Weight: 0.58 kg |
|  OEM103/DG-2 | OEM103 adopts step motor with bracket and shock absorption parts which reduce the vibration and make it easy to load pump head. It has been widely used in flow injection analysis instrument and ultraviolet radiation analysis instrument. | Driver need to be designed or selected to drive the step motor. Acceptable pump heads are YZ and DG series pump heads. | Max. Speed: 100 rpm (DG pump heads) 300 rpm (YZ pump heads) Dimensions: 130 120 115 (mm) |
|  OEMWX15-12 | OEMWX15-12 delivers the maximum flow rates of 46 ml/min, small size; it is suitable to be built in customers equipment. | OEMWX15-12 adopts DC gear motor and #16 silicon tubing, user can control flow rates by adjusting voltage of motor. | Speed: 0-100 rpm, reversible Power supply: DC24V or DC12V Dimension: 70 60 50 (mm) Weight: 0.2kg |

Note: Operating condition is temperature 0 to 40℃, relative humidity <80%.



Features of peristaltic pump tubing

- Good flexibility. Spring back after pressed radially
- Good wear abrasion resistance
- A certain extent pressure bearing capabilities
- Lower gas permeability
- Low absorption, good temperature resistance, not easy to aging, not swelling, anti-corrosion, fewer extractable

Tubing parameters

Inner diameter and wall thickness are the main parameters of the tubing. Different manufacturers have different notations. Such as: specification codes or inner diameter wall thickness.

Tubing materials

Silicon rubber, rubber, plastic, synthetic material, etc. Different materials have different characteristics and different applications.

Tubing selection factors

■ Chemical Compatibility

When transfer different fluids, the tubing should have good chemical characteristics, which is chemical compatibility. Such as: low absorption, good temperature resistance, not easy to aging, not swelling, anti-corrosion, fewer extractable. Chemical resistance decreases as temperature increases. Chemicals that have no effect on the tubing at room temperature could attack the tubing at elevated temperatures.

Chemical compatibility test method:

When the chemical compatibility of the tubing is uncertain, an immersion test is needed. In an immersion test, a small piece of tubing is weighed, and its diameter and length are measured. The tubing then is immersed in a closed vessel with the chemical in question for a minimum of 48 hours. Afterward, the test piece of tubing is rinsed, dried, weighed and measured, and any changes are recorded. The tubing also should be examined for signs of softening or embrittlements, which indicate the chemical has attacked the tubing.

If the tubing survives the test without discoloration, swelling, cracking, loss of flow or other signs of deterioration, then it is compatible with the fluid.

■ Pressure

Peristaltic pump applications typically have been limited by the pressure capabilities of the tubing.

If the pressure is too high, the tubing could swell, resulting in an improper fit through the pump head, which causes excessive wear and tubing failure.

The factors which affect the pressure are material, the proportion of diameter and wall thickness, etc.

■ Temperature

The working temperature range of a tubing material is another important consideration. Different materials have different temperature range.

■ Dimension

The size of the tubing has a direct effect on the amount of fluid delivered. Well-designed pumps are engineered to work with an optimum tubing size or range of tubing sizes, taking into account the tubing's inner diameter and wall thickness. The inner diameter determines the amount of fluid delivered with each turn of the rotor. The wall thickness affects the tubing's ability to spring back to its original shape after each compression, which has a great influence on the overall life of the tubing.

■ Regulatory Approval

Many tubing materials are designed to meet these various regulatory approvals, including those issued by the United States Pharmacopoeia (USP), European Pharmacopoeia (EP), U.S. Food and Drug Administration (FDA), U.S. Department of Agriculture (USDA) and National Sanitation Foundation (NSF).

■ Tubing Flex Life

Different tubing materials have differing abilities to withstand the repeated squeezing action of the rollers. In general, each tubing size, tubing material, pump head style, and operating speed in combination has its own life characteristics. Service life, or flex life in the pump, is the primary concern in a new application. Pump performance is very consistent in a specific application. Maximize the life of a pump system by selecting a tubing material that offers long flex life, using thicker wall tubing, and/or by operating a larger pump at slower speed.

Tubing

Tubing Types which Longer Company Can Provide

Silicon Tubing

■ Platinum-Cured Silicone Tubing

Slightly clarity, smooth surface, low protein binding levels, fewer potential leachable
Ideal for pharmaceutical and biotechnology use, suitable temperature range -51 to 238 °C

■ Peroxide-Cured Silicone Tubing












Slightly clarity, smooth surface, economical, longer tubing life, suitable temperature range -51 to 238°C

SAINT-GOBAIN Tubing












Tygon® , **PharMed® BPT** , **Viton®** , **Fluran®** , **Norprene®** , Please contact Longer Company for detail information

Tubing Specifications:







Medium flow rates tubing

| Tubing Sizes | | 13 [#] | 14 [#] | 19 [#] | 16 [#] | 25 [#] | 17 [#] | 18 [#] | 15 [#] | 24 [#] | 35 [#] | 36 [#] |
|-----------------------------|--------------|---|---|---|---|---|---|---|--|---|---|---|
| Tubing cross sections (1:1) | |  |  |  |  |  |  |  |  |  |  |  |
| Wall thickness (mm) | | 1.6 | | | | | | | 2.4 | | | |
| Inner diameter (mm) | | 0.8 | 1.6 | 2.4 | 3.1 | 4.8 | 6.4 | 7.9 | 4.8 | 6.4 | 7.9 | 9.6 |
| Maximum pressure (Mpa) | Continuous | 0.17 | | | | 0.14 | 0.10 | 0.07 | 0.17 | | 0.14 | |
| | Intermittent | 0.27 | | | | 0.24 | 0.14 | 0.10 | 0.27 | | 0.24 | |

Micro flow rates tubing

| Tubing Sizes | | 0.13 0.86 | 0.5 0.86 | 0.86 0.86 | 1.52 0.86 | 2.06 0.86 | 2.4 0.86 | 2.79 0.86 | 3.17 0.86 | 1 1 | 2 1 | 3 1 |
|-----------------------------|--------------|---|---|---|---|---|---|---|--|---|---|---|
| Tubing cross sections (1:1) | |  |  |  |  |  |  |  |  |  |  |  |
| Wall thickness (mm) | | 0.86 | | | | | | | | | 1.0 | |
| Inner diameter (mm) | | 0.13 | 0.5 | 0.86 | 1.52 | 2.06 | 2.4 | 2.79 | 3.17 | 1.0 | 2.0 | 3.0 |
| Maximum pressure (Mpa) | Continuous | 0.1 | | | | | | | | | | |
| | Intermittent | 0.1 | | | | | | | | | | |

Industrial tubing

| Tubing Sizes | | 73 [#] | 82 [#] | 86 [#] | 90 [#] | 88 [#] | 92 [#] |
|-----------------------------|--------------|---|---|---|--|---|---|
| Tubing cross sections (1:1) | |  |  |  |  |  |  |
| Wall thickness (mm) | | 3.3 | | 6.3 | | 4.8 | |
| Inner diameter (mm) | | 9.6 | 12.7 | 9.6 | 19 | 12.7 | 25.4 |
| Maximum pressure (Mpa) | Continuous | 0.17 | 0.07 | 0.14 | | | |
| | Intermittent | 0.27 | 0.14 | 0.14 | | | |




NAVID KALA DARAN CO.


Tubing Ref. Flow Rates Curve

TUBING


YZ1515x YZII 15




YZII 25




KZ25




YZ2515x



BZ25




YZ35-13 KZ35




Tubing Ref.Flow Rates Curve


DG-1 DG-2 DG-4 DG-6 DG-8 DG-12 (6 rollers)




DG15-28、DG15-48




DG-1 DG-2 DG-4 DG-6 DG-8 DG-12 (10 rollers)




DG15-24



TH15



DMD25



Fitting

Four materials are available for fittings. These fittings match perfectly to the tubing to ensure smooth flow and a proper fit.






High-Density Polyethylene (HDPE): Very good chemical resistance. Temperature range: -53 to 87°C. Sterilize by ethylene oxide only.

Nylon: Good chemical resistance.

Temperature range: -40 to 93°C. Sterilize by ethylene oxide only.

Natural Polypropylene (Natural PP): Very good chemical resistance. Temperature range: -17 to 135°C. Sterilize by ethylene oxide or autoclave.

PVDF: Excellent chemical resistance. Temperature range: -62 to 129°C. Sterilize by ethylene oxide or autoclave.

| Fitting Type | Connector ID | Use with tubing sizes | | |
|---|---------------|-----------------------|-----------------|-------------------------|
| | | Laboratorial Type | Industrial Type | Batch Transferring Type |
|  Straight connectors | 1/16" (1.6mm) | 13# 14# | - | - |
| | 1/8" (3.2mm) | 16# | - | - |
| | 3/16" (4.8mm) | 15# 25# | - | - |
| | 1/4" (6.4mm) | 17# 24# | 26# | - |
| | 3/8" (9.6mm) | 18# 35# 36# | 70# 73# | 86# |
| | 1/2" (12.7mm) | - | 82# 88# | 88# |
| | 5/8" (15.9mm) | - | 89# | - |
| | 3/4" (19.0mm) | - | - | 90# |
| | 1" (25.4mm) | - | - | 92# |
|  "T" type connectors | 1/16" (1.6mm) | 13# 14# | - | - |
| | 1/8" (3.2mm) | 16# | - | - |
| | 3/16" (4.8mm) | 15# 25# | - | - |
| | 1/4" (6.4mm) | 17# 24# | 26# | - |
| | 3/8" (9.6mm) | 18# 35# 36# | 70# 73# | 86# |
| | 1/2" (12.7mm) | - | 82# 88# | 88# |
| | 5/8" (15.9mm) | - | 89# | - |
| | 3/4" (19.0mm) | - | - | 90# |
| | 1" (25.4mm) | - | - | 92# |
|  "Y" type connectors | 1/16" (1.6mm) | 13# 14# | - | - |
| | 1/8" (3.2mm) | 16# | - | - |
| | 3/16" (4.8mm) | 15# 25# | - | - |
| | 1/4" (6.4mm) | 17# 24# | 26# | - |
| | 3/8" (9.6mm) | 18# 35# 36# | 70# 73# | 86# |
| | 1/2" (12.7mm) | - | 82# 88# | 88# |
|  Male threaded connectors | 1/8" - 1/8" | 16# | - | - |
| | 3/16" - 1/8" | 15# 25# | - | - |
| | 1/4" - 1/4" | 17# 24# | 26# | - |
| | 3/8" - 1/4" | 18# 35# 36# | 70# 73# | 86# |
| | 1/2" - 1/2" | - | 82# 88# | 88# |
| | 5/8" - 1/2" | - | 89# | - |
| | 3/4" - 1/2" | - | - | 90# |
|  Reducing connectors | 1/16" - 1/8" | 13# 14# | - | - |
| | 1/8" - 3/32" | 16# | - | - |
| | 3/16" - 1/8" | 15# 25# | - | - |
| | 1/4" - 1/8" | 17# 24# | 26# | - |
| | 3/8" - 1/4" | 18# 35# 36# | 70# 73# | 86# |
| | 1/2" - 3/8" | - | 82# 88# | 88# |
| | 5/8" - 1/2" | - | 89# | - |
| | 3/4" - 1/2" | - | - | 90# |
| | 1" - 1/2" | - | - | 92# |

Accessories

Standard External Control Interface



Applicable Drive Types

BQ50 - 1J BT100 - 1L
BT100 - 2J BT100 - 1F
BT300 - 2J BT300 - 1F
BT600 - 2J WT600 - 2J
LEAD - 2 WT600 - 1F

Five Kinds of Control Modes

0 - 5V analog input
0 - 10V analog input
4 - 20mA analog input
0 - 10kHz Pulse Input
Communication control

Note: Please see above picture for part numbers

Footswitch



| Type | Connector | Part Number | Applicable Products |
|--------|-----------|-------------|--|
| JK - 1 | DB - 9 | 0505201 | Peristaltic Pump |
| JK - 2 | DB - 9 | 0505202 | Dispensing Controller |
| JK - 3 | DB - 15 | 0505203 | Peristaltic Pump |
| JK - 4 | DB - 9 | 0505204 | LSP Series Syringe Pump (the PCB board of the pump needs to be modified) |
| JK - 5 | DB - 15 | 0505205 | TJ Series Syringe Pump |
| JK - 6 | None | 0505203 | TS Series Syringe Pump |

Dispensing Controller



| Type | Part Number | Function and Feature |
|---------|-------------|---|
| FK - 1C | 0505212 | Dispensing controller can be equipped with many peristaltic pumps to realize dispensing automatically |

Footswitch Working Modes

| | |
|------|--|
| JK-1 | 1. Gated: Connect the footswitch to the pump. The pump stops after power up. The pump runs as long as the footswitch is pressed. |
| | 2. Trigger: Connect the footswitch to the pump. The pump stops after power up. Press footswitch, the pump starts running. Press footswitch again, the pump stops. |
| JK-2 | Trigger: Press the footswitch one time, dispensing controller FK-1C starts running according to the set running time. When the set running time is out, FK-1C stops. |
| JK-3 | Gated: The pump runs as long as the footswitch is pressed. |
| JK-4 | 1. Gated: Connect the footswitch to the pump. The pump stops after power up. The pump runs as long as the footswitch is pressed. |
| | 2. Trigger: Connect the footswitch to the pump. The pump stops after power up. Press footswitch, the pump starts running. Press footswitch again, the pump stops. |
| JK-5 | Trigger: Connect the footswitch to the pump. The pump stops after power up. Press footswitch, the pump starts running. Press footswitch again, the pump stops. |
| JK-6 | Trigger: Connect the footswitch to the pump. The pump stops after power up. Press footswitch, the pump starts running. Press footswitch again, the pump stops. |

TJ-1A/L0107-1A SYRINGE PUMP

TJ-2A/L0107-2A



The pumps combine precision, compact size with ease of install and operation. It can hold μL unit standard glass syringe. The features of accurate distance control and broad linear speed range ($7.94\mu\text{m}/\text{min}$ - $79.4\text{mm}/\text{min}$) can meet versatile requirements. Its vertical (horizontal) installation structure makes this pump easily used in micromanipulator, stereotaxic instrument for various biologic research applications.

Functions and Features

Working mode: TJ-1A/L0107-1A: Infusion

TJ-2A/L0107-2A: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous

User-defined glass syringe: Save 4 inner diameters of user-defined glass syringe barrel

Linear force output: Full stroke $>20\text{ N}$

Memory function: Select resume operation or remain stopped when power returns after an interruption

Signal output: Start/Stop output, cw/ccw output (open collector)

Calibration: Acquire accurate volume through calibration

Fast forward & fast reverse: Infusion or filling at the max. Speed

Specifications

Max. Infusion distance: 70 mm

Acceptable glass syringe: 5 - 1000 (μL)

Linear speed: $7.94\mu\text{m}/\text{min}$ - $79.4\text{mm}/\text{min}$

Adjusting resolution: $7.94\mu\text{m}/\text{min}$

Distance resolution: $0.165\mu\text{m}$

Linear force: $>20\text{ N}$

Operating mode: Membrane keypad and rotary encoded switch

Accuracy: $\leq 0.5\%$ error in the condition of $\geq 30\%$ of max. infusion distance

Display: 128×64 graphic LCD

External control: Start/Stop control, fast forward control, fast reverse control

Communication interface: RS485

Power: AC 100 V - 240 V or DC 12 V

Power consumption: $\leq 10\text{ W}$

Operating condition:

Temperature $0 - 40^\circ\text{C}$

Relative humidity $<80\%$

Controller dimensions (L \times W \times H): $170 \times 108 \times 65$ (mm)

Controller weight: 0.8 kg

Drive unit dimensions (L \times W \times H): $180 \times 46 \times 78$ (mm)

Drive unit weight: 0.6 kg

IP rating: IP 21

| Syringe Pump (Part Number) | Glass Syringe Specification (μL) | Barrel Inner Diameter (mm) | Effective Stroke (mm) | Flow Rates (nL/min - $\mu\text{L}/\text{min}$) | Syringe Material | Weight (kg) |
|--|--|-------------------------------|--------------------------|---|---------------------|--|
| TJ - 1A / L0107 - 1A (0503101) (0503001) | 5 | 0.35 | 51.97 | 0.764 - 7.64 | Glass Syringe | Controller 0.8 Drive Unit 0.6 |
| | 10 | 0.50 | 50.93 | 1.559 - 15.59 | | |
| | 25 | 0.80 | 49.74 | 3.989 - 39.89 | | |
| | 50 | 1.10 | 52.61 | 7.544 - 75.44 | | |
| TJ - 2A / L0107 - 2A (0503111) (0503011) | 100 | 1.60 | 49.74 | 15.96 - 159.6 | | |
| | 250 | 2.30 | 60.17 | 32.98 - 329.8 | | |
| | 500 | 3.25 | 60.27 | 65.85 - 658.5 | | |
| | 1000 | 4.61 | 59.91 | 132.5 - 1325 | | |

Note: Flow rates = Linear rate \times Section area of the barrel

TS-1A/L0107-1A SYRINGE PUMP

TS-2A/L0107-2A



The pumps combine precision, compact size, multiple functions with ease of operation. It can hold μL unit standard glass syringe. The features of accurate distance control and broad linear speed range ($7.94\mu\text{m}/\text{min}$ - $79.4\text{mm}/\text{min}$) can meet versatile requirements. The drive unit is separate, easy to installation and combination. Its vertical (horizontal) installation structure makes this pump easily used in micromanipulator, stereotaxic instrument for various biologic research applications.

TS-2A/L0107-2A
Infusion/Withdrawal
Lock devices
See **View**

Specifications

Max. infusion distance: 70 mm
Acceptable glass syringe: $5\mu\text{L}$ - $1000\mu\text{L}$
Linear speed: $7.94\mu\text{m}/\text{min}$ - $79.4\text{mm}/\text{min}$
Adjusting resolution: $7.94\mu\text{m}/\text{min}$
Distance resolution: $0.165\mu\text{m}$
Linear force: $>20\text{N}$
Operating mode: Membrane keypad and rotary encoded switch
Accuracy: $\leq 0.5\%$ error in the condition of $\geq 30\%$ of max. infusion distance
Display: 128×64 graphic LCD
External control: Start/Stop control, fast forward control, fast reverse control
Communication interface: RS485
Power: AC 100 V - 240 V or DC 12 V
Power consumption: $<40\text{W}$
Operating condition: Temperature $0 - 40^\circ\text{C}$
 Relative humidity $<80\%$
Controller dimensions (L \times W \times H):
 $170 \times 108 \times 65$ (mm)
Controller weight: 0.9 kg
Drive unit dimensions (L \times W \times H):
 $180 \times 46 \times 78$ (mm)
Drive unit weight: 0.6 kg
IP rating: IP 21

Functions and Features

Parameters setting: The parameters of each channel can be different.
Running control: Each drive unit can be controlled separately; Or four drive units can be controlled to run simultaneously or run at different time.
Channel copy: All drive units can run according to the parameters of one drive units of them
Delaying startup: Delaying startup time of each channel can be controlled separately
Memory function: Select resume operator or remain stopped when power returns after an interruption
Block protection: When one drive unit stops accidentally, the system will warm and stop
Working mode: TS-1A/L0107-1A: Infusion
 TS-2A/L0107-2A: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous
External control: Start/stop input control signal which is pulse mode to swith the states of start and stop
 Each channel has two ways OC gate output to indicate the start/stop and direction of the channel
Communication: Realize computer control through RS485 communication interface

Main Functions for Each Channel

Syringe selection: The syringe can be selected in the manufacturer table which includes manufacturer, material and size
User-defined glass syringe: Save 4 inner diameters of user-defined glass syringe barrel
Parameters setting: Set dispensing volume, infusion time, pause time and copy number
Display mode selection: Different parameters (volume, flow rate, linear speed) can be selected in the main display interface
Fast forward & fast reverse: Infusion or withdrawal at the max. speed
Calibration: Acquire accurate volume through calibration

| Syringe Pump (Part Number) | Glass Syringe Specification (μL) | Barrel Inner Diameter (mm) | Effective Stroke (mm) | Flow Rates (nL/min - $\mu\text{L}/\text{min}$) | Syringe Material | Weight (kg) |
|--|--|-------------------------------|--------------------------|---|---------------------|--|
| TS - 1A / L0107 - 1A (0503151) (0503001) | 5 | 0.35 | 51.97 | 0.764 - 7.64 | Glass Syringe | Controller 0.8 Drive Unit 0.6 |
| | 10 | 0.50 | 50.93 | 1.559 - 15.59 | | |
| | 25 | 0.80 | 49.74 | 3.989 - 39.89 | | |
| | 50 | 1.10 | 52.61 | 7.544 - 75.44 | | |
| TS - 2A / L0107 - 2A (0503161) (0503011) | 100 | 1.60 | 49.74 | 15.96 - 159.6 | | |
| | 250 | 2.30 | 60.17 | 32.98 - 329.8 | | |
| | 500 | 3.25 | 60.27 | 65.85 - 658.5 | | |
| | 1000 | 4.61 | 59.91 | 132.5 - 1325 | | |

Note: Flow rates=Linear rate \times Section area of the barrel

TJ-3A/W0109-1B Syringe Pump



This pump is a single channel micro syringe pump which has infusion/withdrawal mode. It combines precision, compact size, multiple functions with ease of operation. It can hold milliliter unit standard syringe. The features of accurate distance control and broad linear speed range (7.94 $\mu\text{m}/\text{min}$ - 79.4 mm/min) can meet versatile requirements. The drive unit is independent. It is easy to combination and installation. It can meet different operation requirements and is suitable for various research fields.

Functions and Features

Working mode: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous

User-defined syringe: Save 4 inner diameters of user-defined syringe barrel

Linear force output: Whole course $>90\text{ N}$

Memory function: Select resume operation or remain stopped when power returns after an interruption

Signal output: Start/Stop output, cw/ccw output (open collector)

Calibration: Acquire accurate volume through calibration

Fast forward & fast reverse: Infusion or withdrawal at the max. speed

Specifications

Max. infusion distance: 90 mm

Acceptable syringe: 5 μL - 60 mL

Linear speed: 7.94 $\mu\text{m}/\text{min}$ - 79.4 mm/min

Adjusting resolution: 7.94 $\mu\text{m}/\text{min}$

Distance resolution: 0.165 μm

Linear force: $>90\text{ N}$

Operating mode: zMembrane keypad and rotary encoded switch

Accuracy: $\leq \pm 0.5\%$ error in the condition of $\geq 30\%$ of max. infusion distance

Display: 128 \times 64 graphic LCD

External control: Start/stop control, fast forward control, fast reverse control

Communication interface: RS485

Power: AC 90 V - 240 V or DC 12 V

Power consumption: $\leq 10\text{ W}$

Operating condition: Temperature 0 - 40 $^{\circ}\text{C}$

Relative humidity $< 80\%$

Controller dimensions (L \times W \times H): 170 \times 108 \times 65 (mm)

Controller weight: 0.8 kg

Drive unit dimensions (L \times W \times H): 245 \times 100 \times 95 (mm)

Drive unit weight: 1.3 kg

| Syringe Pump (Part Number) | Glass Syringe Specification (mL) | Barrel Inner Diameter (mm) | Effective Stroke (mm) | Flow Rates ($\mu\text{L}/\text{min}$ - mL/min) | Syringe Material | Weight (kg) |
|--|-------------------------------------|-------------------------------|--------------------------|---|----------------------------------|--|
| TJ - 3A / W0109 - 1B (0503121) (0503021) | 1 | 4.7 | 57.0 | 0.138 - 1.38 | Class Syringe Plastic Syringe | Controller 0.8 Drive Unit 1.3 |
| | 2 | 8.98 | 31.6 | 0.503 - 5.03 | | |
| | 5 | 12.25 | 42.5 | 0.936 - 9.36 | | |
| | 10 | 14.9 | 58.2 | 1.384 - 13.84 | | |
| | 20 | 19.05 | 70.0 | 2.262 - 22.62 | | |
| | 30 | 22.05 | 78.0 | 3.031 - 30.31 | | |
| | 60 | 29.15 | 89.9 | 5.297 - 52.97 | | |

Note: Flow rates = Linear rate \times Section area of the barrel

TS-1B/W0109-1B Syringe Pump

This pump is a four-channel micro syringe pump which has infusion/withdrawal mode. It combines precision, compact size, multiple functions with ease of operation. It can hold milliliter unit standard syringe. The features of accurate distance control and broad linear speed range (7.94 $\mu\text{m}/\text{min}$ - 79.4 mm/min) can meet versatile requirements. The drive unit is independent. It is easy to combination and installation. It can perform complex operation and is suitable for various research fields.



Specifications

Max. infusion distance: 90 mm
Acceptable syringe: 5 μL - 60 mL
Linear speed: 7.94 $\mu\text{m}/\text{min}$ - 79.4 mm/min
Adjusting resolution: 7.94 $\mu\text{m}/\text{min}$
Distance resolution: 0.165 μm
Linear force: $>90\text{ N}$
Operating mode: Membrane keypad and rotary encoded switch
Accuracy: $\leq 0.5\%$ error in the condition of $\geq 30\%$ of max. infusion distance
Display: 128 \times 64 graphic LCD
External control: Start/stop control
Communication interface: RS485
Power: AC 90 V - 240 V or DC 12 V
Power consumption: $\leq 40\text{ W}$
Operating condition:
 Temperature 0 - 40 $^{\circ}\text{C}$
 Relative humidity $<80\%$
Controller dimensions (L \times W \times H):
 235 \times 178 \times 74 (mm)
Controller weight: 0.9 kg
Drive unit dimensions (L \times W \times H):
 245 \times 100 \times 95 (mm)
Drive unit weight: 1.3 kg

Functions and Features

Parameters setting: The parameters of each channel can be different.
Running control: Each drive unit can be controlled separately; Or four drive units can be controlled to run simultaneously or run at different time.
Channel copy: All drive units can run according to the parameters of one drive units of them
Delaying startup: Delaying startup time of each channel can be controlled separately
Memory function: Select resume operation or remain stopped when power returns after an interruption
Block protection: When one drive unit stops accidentally, the system will warn and stop
Working mode: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous
External control: Accept pulse signal to control the start/stop of the pump. Each channel has two ways OC gate output to indicate the start/stop and direction of the channel
Communication: Realize computer control through RS485 communication interface

Main Functions for Each Channel

Syringe selection: The syringe can be selected in the manufacturer table which includes manufacturer, material and size
User-defined syringe: Save 4 inner diameters of user-defined glass syringe barrel
Working mode: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous
Parameters setting: Set dispensing volume, infusion time, withdrawal time, pause time and copy number
Display mode selection: Different parameters (volume, flow rates, linear speed) can be selected in the main display interface
Fast forward & fast reverse: Infusion and withdrawal at the max. speed
Calibration: Acquire accurate volume through calibration

| Syringe Pump (Part Number) | Glass Syringe Specification (mL) | Barrel Inner Diameter (mm) | Effective Stroke (mm) | Flow Rates ($\mu\text{L}/\text{min}$ - mL/min) | Syringe Material | Weight (kg) |
|--|-------------------------------------|-------------------------------|--------------------------|---|----------------------------------|--|
| TS - 1B / W0109 - 1B (0503171) (0503021) | 1 | 4.7 | 57.0 | 0.138 - 1.38 | Class Syringe Plastic Syringe | Controller 0.8 Drive Unit 1.3 |
| | 2 | 8.98 | 31.6 | 0.503 - 5.03 | | |
| | 5 | 12.25 | 42.5 | 0.936 - 9.36 | | |
| | 10 | 14.9 | 58.2 | 1.384 - 13.84 | | |
| | 20 | 19.05 | 70.0 | 2.262 - 22.62 | | |
| | 30 | 22.05 | 78.0 | 3.031 - 30.31 | | |
| | 60 | 29.15 | 89.9 | 5.297 - 52.97 | | |

Note: Flow rates = Linear rate \times Section area of the barrel

LSP01-1A/2A, LSP04-1A Syringe Pump

LSP01-1A/2A



LSP04-1A



Functions and Features

Syringe selection: The syringe can be selected in the manufacturer table which includes manufacturer, material and size or input the inner diameter of the syringe barrel directly

Easy to operate: Combining big screen LCD display with rotary encoded switch and membrane keypad makes the operation simple and prompt

Working mode: Infusion

Memory function:

1. The parameters are saved in EEPROM. The parameters don't need to be reset when power returns after an interruption

2. In flow rates mode, the pump remains running or stop according to the setting parameters when power returns after an interruption

Protection function: The pump will stall and give an alarm when the drive structure of the pump is blocked

Communication function: Realize computer control through RS485 communication interface

External control function: Input/output control

Calibration function: Acquire accurate volume through calibration

Syringe protection: Adjust syringe rest to prevent syringe from damaging

Infusion only syringe pumps. The acceptable syringes for LSP01-1A/2A are from 10 μ L to 60 mL and for LSP04-1A is from 10 μ L to 10 mL. Suitable for high accuracy and small flow rates liquid transferring.

Specifications Comparison Table

| Syringe Pump | LSP01 - 1A | LSP01 - 2A | LSP04 - 1A |
|-------------------------------|--|-------------------------------|------------------------------------|
| Max. No. of syringes | 1 | 1 | 4 |
| Infusion volume per microstep | 0.13 μ L (60 mL Syringe) | 0.026 μ L (60 mL Syringe) | 0.026 μ L (10mLSyringe) |
| Syringe size | 10 μ L - 60mL | 10 μ L - 60mL | 10 μ L - 10mL |
| Flow rates | 0.831nL/min - 54.155mL/min | 0.166nL/min - 10.83mL/min | 0.831nL/min - 21.675mL/min |
| Advance per microstep | 1/16step: 0.156 μ m | 1/16step: 0.03125 μ m | 1/16step: 0.156 μ m |
| Max. linear rate | 65mm/min | 13mm/min | 130mm/min |
| Min. linear rate | 5 μ m/min | 1 μ m/min | 5 μ m/min |
| Max. step rate | 6933 (1/16step) / sec | | 13867 (1/16step) /sec |
| Min. step rate | 16 (1/16 step) /30 sec | | |
| Working mode | Infusion | | |
| Linear force | > 9kgf | | |
| Accuracy | $\leq \pm 0.5\%$ error in the condition of > 30% of max. infusion distance | | |
| Operating mode | Rotary encoded switch and membrane keypad | | |
| Display | 128 x 64 graphic LCD | | |
| Power | AC 100 - 240 V | | |
| Operating condition | Temperature 5 $^{\circ}$ C - 40 $^{\circ}$ C Relative humidity < 80 % | | |
| Dimensions | 280 \times 210 \times 140 (mm) | | 280 \times 250 \times 140 (mm) |
| Weight | 3.6 kg | | 4.5kg |

| Syringe Pump | Syringe | I.D. (mm) | Flow Rates (μ L/min - mL/min) | Linear Rate | Part Number |
|--------------|------------|-----------|---------------------------------------|---------------------------|-------------|
| LSP01 - 1A | 50 μ L | 1.03 | 0.004 - 0.054 | 5 μ m/min - 65mm/min | 0503401 |
| | 10mL | 14.57 | 0.834 - 10.837 | | |
| | 60mL | 32.57 | 4.166 - 54.155 | | |
| LSP01 - 2A | 50 μ L | 1.03 | 0.001 - 0.011 | 1 μ m/min - 13mm/min | 0503411 |
| | 10mL | 14.57 | 0.167 - 2.167 | | |
| | 60mL | 32.57 | 0.833 - 10.831 | | |
| LSP04 - 1A | 50 μ L | 1.03 | 0.004 - 0.108 | 5 μ m/min - 130mm/min | 0503451 |
| | 1mL | 4.61 | 0.083 - 2.170 | | |
| | 10mL | 14.57 | 0.834 - 21.675 | | |

Note: Flow rates=Linear rate \times Section area of the barrel

LSP02/10-1B Syringe Pump

These pumps are multi-channel syringe pump which have infusion/withdrawal mode. The acceptable syringes for LSP02-1B are from 10 μ L to 140mL and the acceptable syringes for LSP10-1B are from 10 μ L to 10mL. Suitable for high accuracy and small flow rates liquid transferring.

Functions and Features

Syringe selection: The syringe can be selected in the manufacturer table which includes manufacturer, material and size or input the inner diameter of the syringe barrel directly

Easy to operate: Combining big screen LCD display with rotary encoded switch and membrane keypad makes the operation simple and prompt.

Working mode: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous

Memory function:

1. The parameters are saved in EEPROM. The parameters don't need to be reset when power returns after an interruption.
2. In flow rates mode, the pump remains running or stop according to the setting parameters when power returns after an interruption.

Protection function: The pump will stall and give an alarm when the drive structure of the pump is blocked

Communication function: Realize computer control through RS485 communication interface

External control function: Input/output control

Calibration function: Acquire accurate volume through calibration

Syringe protection: Adjust syringe rest to prevent syringe from damaging



LSP02-1B



LSP10-1B

Specifications Comparison Table

| Syringe Pump | LSP02 - 1B | LSP10 - 1B |
|-------------------------------|--|------------------------------------|
| Max. No. of syringes | 2 | 10 |
| Infusion volume per microstep | 0.13 μ L (60 mL BD Syringe) | 0.0257 μ L (10mL BD Syringe) |
| Syringe size | 10 μ L - 140mL | 10 μ L - 10mL |
| Flow rates | 0.831nL/min - 150.5mL/min | 0.831nL/min - 21.675mL/min |
| Advance per microstep | 1/16step: 0.156 μ m | |
| Max. step rate | 13867 (1/16step) /sec | |
| Min. step rate | 16 (1/16 step) /30sec | |
| Max. linear rate | 130mm/min | |
| Min. linear rate | 5 μ m/min | |
| Working mode | Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous | |
| Linear force | >18kgf | |
| Accuracy | $\leq \pm 0.5\%$ error in the condition of >30% of max. infusion distance | |
| Operating mode | Rotary encoded switch and membrane keypad | |
| Display | 128 \times 64 graphic LCD | |
| Power | AC 100 - 240 V | |
| Operating condition | Temperature 5 $^{\circ}$ C-40 $^{\circ}$ C Relative humidity: < 80% | |
| Dimensions | 280 \times 250 \times 140 (mm) | 280 \times 330 \times 140 (mm) |
| Weight | 4.3 kg | 5.3 kg |

| Syringe | I.D. (mm) | LSP02 - 1B | | | LSP10 - 1B | | |
|------------|-----------|--|--|-------------|---------------------------------------|--|-------------|
| | | Ref. Flow Rates (μ L/min - mL/min) | Linear Rate (μ m/min - mm/min) | Part Number | Flow Rates (μ L/min - mL/min) | Linear Rate (μ m/min - mm/min) | Part Number |
| 50 μ L | 1.03 | 0.004 - 0.108 | 5 - 130 | 0503441 | 0.004 - 0.108 | 5-130 | 0503461 |
| 1mL | 4.61 | 0.083 - 2.170 | | | 0.083 - 2.170 | | |
| 2.5mL | 7.28 | 0.208 - 5.411 | | | 0.208 - 5.411 | | |
| 5mL | 10.30 | 0.417 - 10.832 | | | 0.417 - 10.832 | | |
| 10mL | 14.57 | 0.834 - 21.675 | | | 0.834 - 21.675 | | |
| 25mL | 23.03 | 2.083 - 54.153 | | | | | |
| 50mL | 32.57 | 4.166 - 108.310 | | | | | |
| 60mL | 32.57 | 4.166 - 108.310 | | | | | |

Note: Flow rates=Linear rate \times Section area of the barrel

LSP01-1C Syringe Pump



This pump is a two-syringe push-pull syringe pump. The acceptable syringes are from 10 μ L to 10mL. Suitable for high accuracy and small flow rate liquid continuously transferring.

Functions and Features

Syringe selection: The syringe can be selected in the manufacturer table which includes manufacturer, material and size or input the inner diameter of the syringe barrel directly

Easy to operate: Combining big screen LCD display with rotary encoder switch and membrane keypad makes the operation simple and prompt

Working mode: Push-pull

Memory function:

1. The parameters are saved in EEPROM. The parameters don't need to be reset when power returns after an interruption
2. In flow rate mode, the pump remains running or stop according to the setting parameters when power returns after an interruption

Protection function: The pump will stall and give an alarm when the drive structure of the pump is blocked

Communication function: Realize computer control through RS485 communication interface

External control function: Input/output control

Calibration function: Acquire accurate volume through calibration

Syringe protection: Adjust syringe rest to prevent syringe from damaging

Specifications

| | |
|-------------------------------|--|
| Syringe size | 10 μ L - 10mL |
| Linear force | 9kg |
| Advance per microstep | 0.156 μ m (1/16step) |
| Infusion volume per Microstep | 0.026 μ L (10mL syringe 1/16step) |
| Max. step rate | 6933 (1/16step)/sec |
| Min. step rate | 16 (1/16step)/30sec |
| Max. linear rate | 65mm/min |
| Min. linear rate | 5 μ m/min |
| Flow rates | 0.831nL/min - 10.84mL/min |
| Accuracy | \leq 0.5% error in the condition of \geq 30% of max. Infusion distance |
| Operating mode | Rotary coded switch and membrane keypad |
| Display | 128 \times 64 graphic LCD |
| Power | AC100-240V |
| Operating condition | Temperature 5 $^{\circ}$ C - 40 $^{\circ}$ C Relative humidity <80% |
| Dimensions | 280 \times 220 \times 140 (mm) (L \times W \times H) |
| Weight | 3.6kg |

| Syringe Pump | Part Number | Syringe | Inner Diameter (mm) | Flow Rates | Linear Rate | Weight (kg) |
|--------------|-------------|------------|---------------------|--------------------------------------|--------------------------|-------------|
| LSP01 - 1C | 0503421 | 10 μ L | 0.46 | 0.049 μ L/hr - 10.80 μ L/min | 5 μ m/min - 65mm/min | 3.6 |
| | | 1mL | 4.61 | 5.007 μ L/hr - 1085 μ L/min | | |
| | | 2.5mL | 7.28 | 12.49 μ L/hr - 2706 μ L/min | | |
| | | 5mL | 10.30 | 25.00 μ L/hr - 5415 μ L/min | | |
| | | 10mL | 14.57 | 50.02 μ L/hr - 10833 μ L/min | | |

Note: Flow rates=Linear rate \times Section area of the barrel

LSP01-1BH Syringe Pump



This pump is a single channel high pressure syringe pump. The acceptable syringes are 2.5mL-100mL stainless steel syringes. Because of the extra linear force this pump is suitable for transferring viscous fluids or delivering the fluid to reactors in chemical applications.

Functions and Features

Syringe selection: The syringes can be selected in the manufacturer table which includes manufacturer, material and size or input the inner diameter of the syringe barrel directly

Easy to operate: Combining big screen LCD display with rotary encoder switch and membrane keypad makes the operation simple and prompt.

Working mode: Infusion, withdrawal, infusion/withdrawal, withdrawal/infusion, continuous

Memory function:

1. The parameters are saved in EEPROM . The parameters don't need to be reset when power returns after an interruption
2. In flow rate mode, the pump remains running or stop according to the setting parameters when power returns after an interruption

Protection function: The pump will stall and give an alarm when the drive structure of the pump is blocked

Communication function: Realize computer control through RS485 communication interface

External control function: Input/output control

Calibration function: Acquire accurate volume through calibration

Syringe protection: Adjust syringe rest to prevent syringe from damaging

Specifications


| | |
|-------------------------------|--|
| Syringe size | 2.5 - 100 (mL) |
| Linear force | 45kg |
| Advance per microstep | 0.156 μ m (1/16step) |
| Infusion volume per microstep | 0.149 μ L (100mL syringe 1/16step) |
| Max. Step rate | 13867 (1/16step) /sec |
| Min. step rate | 16 (1/16step) /30sec |
| Max. linear rate | 130mm/min |
| Min. Linear rate | 5 μ m/min |
| Flow rates | 0.09 μ L/min - 124.36mL/min |
| Accuracy | \leq 0.5% error in the condition of \geq 30% of max. Infusion distance |
| Operating mode | Rotary coded switch and membrane keypad |
| Display | 128 64 graphic LCD |
| Power | AC100 - 240V |
| Operating condition | Temperature 5 $^{\circ}$ C - 40 $^{\circ}$ C Relative humidity <80% |
| Dimensions | 280 250 140 (mm) (L \times W \times H) |
| Weight | 6.3kg |

| Syringe Pump | Part Number | Inner Diameter (mm) | Syringe (mL) | Flow Rates | Outlet Pressure | Linear Rate | Weight (kg) |
|--------------|-------------|---------------------|--------------|----------------------------------|-----------------|---------------------------|-------------|
| LSP01 - 1BH | 0503431 | 4.79 | 2.5 | 5.406 μ L/hr - 140.5mL/hr | >19.48MPa | 5 μ m/min - 130mm/min | 6.3 |
| | | 9.53 | 8 | 21.399 μ L/hr - 556.3mL/hr | >5.05MPa | | |
| | | 19.13 | 20 | 86.226 μ L/hr - 2241.8mL/hr | >1.25MPa | | |
| | | 28.6 | 50 | 192.727 μ L/hr - 5010.9mL/hr | >0.56MPa | | |
| | | 34.9 | 100 | 286.986 μ L/hr - 7461.6mL/hr | >0.37MPa | | |

Note: Flow rates=Linear rate \times Section area of the barrel

WT3000-1JA Micro Gear Pump

Flow Rates



Applicable Pump Heads

MG204
MG209
MG213



WT3000-1JA micro gear pump adopts brushless DC motor which features high efficiency and free maintenance. It can deliver flow rates from 85.7 to 2571.4 ml/min with the features of small pulse and big pressure.

Specifications

| | |
|--|--|
| Flow rates | 85.7-2571.4 (mL/min) |
| Speed | 300-3000 (rpm) |
| Speed precision | 1.0% |
| Speed resolution | 1.0 rpm |
| Viscosity of liquid | ≤200 cSt |
| External input control | Speed control (4-20mA, 0.5-5V, 1.0-10V, 1.0-10kHz corresponding to 300-3000 rpm), start/stop control |
| Diameter of particle in liquid | ≤10μm |
| Liquid temperature range | -45-50℃ with PTFE gear, -45-120℃ with PEEK gear |
| Communication interface | RS485 |
| Max. outlet pressure (testing medium is water) | MG204: 0.8MPa MG209: 0.8MPa MG213: 0.3MPa |
| Power supply | AC 176-264V/90-130V, 50/60 Hz |
| Power consumption | ≤50W |
| Operating condition | Temperature 0 to 40℃ Relative humidity <90% (non-condensing) |
| Drive dimensions | 232×142×149 (mm) (L×W×H) |
| Weight | 2.83 kg |
| IP rating | IP31 |


Functions and Features

- Acceptable pump head: MG204, MG209, MG213
- Operating mode: Membrane keypad
- Prime function for fast filling and emptying
- Display: 4 digits LED displays the speed of the pump, two LED lights indicate the running status and external control status
- External input function: Control the start/stop and the speed of the pump
- Footswitch: Control the start/stop of the pump
- Communication function: Realize computer control through RS485 interface
- Memory function: Store the running parameters automatically
- Cooling mode: Heat-emitting fan
- Running direction: Clockwise

| Drive | Part Number | Pump Head | Speed (rpm) | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|-----------|-------------|---------------------|-------------|
| WT3000-1JA | 0502701 | MG204 | 300-3000 | 85.7-857.1 | 3.25 |
| | | MG209 | | 171.4-1714.3 | 3.26 |
| | | MG213 | | 257.1-2571.4 | 3.28 |

WT3000-1JB Micro Gear Pump

Flow Rates




Applicable Pump Heads

MS204
MS209
MS213



Pump Head Structure Drawing



1.Pump Head 2.Outer Magnet 3.Adapter Kit

WT3000-1JB adopts brushless DC motor. It has features of high efficiency and free maintenance. It accepts multiple pump heads and delivers flow rates from 85.7 to 2571.4 mL/min. The compact and precise gear pump heads are noiseless, pulsation-free and suitable to transfer viscous, high temperature fluids or transfer the fluids in pressure condition.

Specifications

| | |
|---|--|
| Flow rates | 85.7 - 2571.4 (mL/min) |
| Speed | 300 - 3000 (rpm) |
| Speed accuracy | 1.0% |
| Speed resolution | 1.0 rpm |
| Fluids viscosity | ≤200 cSt |
| External control | Start/stop control, speed control (0.5 - 5V, 1 - 10 V, 4 - 20 mA, 1 - 10 kHz optional) |
| Communication interface | RS485 |
| Max. outlet pressure (testing medium is water) | MS204: 1.4MPa MS209: 0.9MPa MS213: 0.8MPa |
| Power supply | AC 220/110V ±20%, 50Hz/60Hz |
| Power consumption | ≤150 W |
| Operating condition | Temperature 0 to 40℃ Relative humidity <90% (non-condensing) |
| Drive weights | 5.1 kg |
| Drive dimensions | 290×207×180 (mm) (L×W×H) |
| IP rating | IP 31 |

Functions and Features


- Acceptable pump heads: MS204, MS209, MS213
- Suitable to transfer viscous, high temperature fluids or transfer the fluids in pressure condition
- Compact size, noiseless, pulsation-free
- 4 digits LED displays the speed of the pump
- Combine the knob and rotary encoder switch to make the operation easy
- Prime function for fast filling and emptying
- The start/stop and speed of the pump can be controlled through external control port
- Store the running parameters automatically
- Decoupling will occur when the pump load exceeds the maximum coupling torque provided by the alignment of the two magnets

| Drive | Part Number | Pump Head | Speed (rpm) | Flow Rates (mL/min) | Weight (kg) |
|--------------|-------------|-----------|-------------|---------------------|-------------|
| WT3000 - 1JB | 0502702 | MS204 | 300 - 3000 | 85.7 - 857.1 | 5.39 |
| | | MS209 | | 171.4 - 1714.3 | 5.40 |
| | | MS213 | | 257.1 - 2571.4 | 5.42 |

WT3000-1FA Micro Gear Pump



Flow Rates



Applicable Pump Heads

MG204
MG209
MG213



WT3000-1FA micro gear pump has two work modes of flow rate and dispensing to deliver flow rates from 85.7 to 2571.4 ml/min and dispense liquid volume from 0.1 ml to 999 liter. It adopts 128×32 LCD display all the running parameters. Stainless pump head is compact and precise with the features of low noise and small pulse. It is suitable to deliver high pressure and high temperature liquid.

Specifications

| | |
|---|---|
| Flow rates | 85.7-2571.4 (mL/min) |
| Dispensing volume | 0.1mL-999L |
| Copy No. | 0-9999 , '0' means running continuously |
| Pause time | 1s-999h, min. adjustment resolution: 0.1s |
| Flow rate calibration time | 30-1800 second, adjustment resolution: 1s |
| Viscosity of liquid | ≤200 cSt |
| External input control | Flow rate control (4-20mA, 0.5-5V, 1.0-10V, 1.0-10kHz) , start/stop control |
| External output control | 1.25-12.5kHz frequency output corresponding to 300-3000 rpm; start/stop status output |
| Diameter of particle in liquid | ≤10μm |
| Max. outlet pressure (testing medium: water) | MG204: 0.8MPa MG209: 0.8MPa MG213: 0.3MPa |
| Liquid temperature range | -45-50℃ with PTFE gear, -45-120℃ with PEEK gear |
| Communication interface | RS485 |
| Power supply | AC 176-264V/90-130V, 50/60 Hz |
| Power consumption | ≤50W |
| Operating condition | Temperature 0 to 40℃ Relative humidity <90% (non-condensing) |
| Drive dimensions | 232×142×149 (mm) (L×W×H) |
| Weight | 2.83 kg |
| IP rating | IP31 |


Functions and Features

- Acceptable pump heads: MG204, MG209, MG213
- Flow rate function: Deliver the liquid according to the set flow rates
- Dispensing function: Dispensing volume, copy number and pause time can be set
- Calibration function: The flow rates and the dispensing volume can be calibrated to increase the accuracy
- Display: 128 32 LCD display all the running parameters
- Operating mode: Membrane keypad
- Prime function for fast filling and emptying
- External control input: Control the flow rates and start/stop of the pump
- External control output: Output start/stop and speed signal
- Footswitch: Control the start/stop of the pump
- Communication function: Realize computer control through RS485 interface
- Memory function: Store the running parameters automatically
- Cooling mode: Heat-emitting fan
- Running direction: Clockwise

| Drive | Part Number | Pump Head | Speed (rpm) | Flow Rates (mL/min) | Weight (kg) |
|------------|-------------|-----------|-------------|---------------------|-------------|
| WT3000-1FA | 0502711 | MG204 | 300-3000 | 85.7 - 857.1 | 3.25 |
| | | MG209 | | 171.4 - 1714.3 | 3.26 |
| | | MG213 | | 257.1 - 2571.4 | 3.28 |

WT3000-1FB Micro Gear Pump

Flow Rates




Applicable Pump Heads

MS204
MS209
MS213



Pump Head Installation Drawing



WT3000-1FB gear pump has two work modes of flow mode and dispensing mode. It delivers flow rates from 85.7 to 2571.4 mL/min and the dispensing volume is from 0.1 mL to 999L. 128×32 LCD displays the parameters and information. The compact and precise gear pump heads are noiseless, pulsation-free and suitable to transfer high temperature and high pressure fluids.

Specifications

| | |
|---|---|
| Flow rates | 85.7 - 2571.4 (mL/min) |
| Dispensing volume | 0.1 mL - 999 L |
| Copy number | 0 - 9999 ("0" means running continuously) |
| Pause time | 1s - 999 h |
| Fluid viscosity | ≤200 cSt |
| External control input | Start/stop control, flow rates control (0.5 - 5 V, 1 - 10V, 4 - 20 mA, 1 - 10 kHz optional) |
| External control output | Pulse output 0 - 12.5 kHz (corresponding to 0-3000 rpm) , start/stop state output |
| Communication interface | RS485 |
| Max. outlet pressure (testing medium is water) | MS204: 1.4MPa MS209: 0.9MPa MS213: 0.8MPa |
| Power supply | AC 220/110V ±20%, 50Hz/60Hz |
| Power consumption | ≤150 W |
| Operating condition | Temperature 0 to 40℃ Relative humidity <90% (non-condensing) |
| Drive weights | 5.1 kg |
| Drive Dimensions | 290×207×180 (mm) (L×W×H) |
| IP rating | IP 31 |

Functions and Features

- Acceptable pump heads: MS204, MS209, MS213
- Suitable to transfer viscous, high temperature fluids or transfer the fluids in pressure condition
- Compact size, noiseless, pulsation-free
- Dispensing volume, copy number and pause time can be set
- Flow rates and dispensing volume can be calibrated for higher accuracy
- Combine the membrane keypad and rotary encoder switch to make the operation easy
- Prime function for fast filling and emptying
- The pump has external control input and output ports
- Store the running parameters automatically
- Decoupling will occur when the pump load exceeds the maximum coupling torque provided by the alignment of the two magnets


| Drive | Part Number | Pump Head | Speed (rpm) | Flow Rates (mL/min) | Weight (kg) |
|--------------|-------------|-----------|-------------|---------------------|-------------|
| WT3000 - 1FB | 0502712 | MS204 | 300-3000 | 85.7 - 857.1 | 5.39 |
| | | MS209 | | 171.4 - 1714.3 | 5.40 |
| | | MS213 | | 257.1 - 2571.4 | 5.42 |

SMD02-1 Diluters & Dispensers System



SMD02-1 Dispenser & Dilutor handles liquid sample precisely. It consists of software, two pieces control units and probe. It has functions of semi-automatically dispensing and diluting. SMD02-1 features accurate control, easy and reliable operation to match the requirements to prepare liquid sample for analytical laboratory. SMD02-1 completely replaces the traditional manual operation of dispensing and diluting and widely used in blood analysis, emulsion analysis, chemical analysis, food analysis, urine analysis.

Accessories




SMD01-1 Diluters & Dispensers System

SMD01-1 Dispenser & Dilutor has one piece control unit. Other features and functions are the same with SMD02-1 except dilution proportion range.



Accessories



Diluters & Dispensers System

Functions and Features

- User-friendly interface. Easy to operate.
- The precision of handling volume is $\leq \pm 0.5\%$, and the accuracy of handling volume is $\leq 99\%$ that avoid the artificial error between different operators.
- Save the preparation time of testing and the consumption of solvent.
- Five dispensing modes are available.
- Communication: Bluetooth wireless and RS232.

Work Modes

Dispensing mode:

1. Withdraw liquid to syringe according to programmed speed and volume.
2. Dispense liquid from syringe according to programmed speed and volume.

Diluting mode:

1. Withdrawing liquid from solvent bottle according to programmed parameters.
2. Withdraw liquid from sample bottle to probe.
3. Dispense sample and solvent to bottle.

User define mode:

According to customer actual need, there are four operation modes for each syringe (inlet withdrawing, inlet dispensing, outlet withdrawing, outlet dispensing) to complete complicate work.

Technical Specifications

| Type | SMD02-1 | | SMD01-1 |
|----------------|--------------------------|---|----------|
| Specifications | Acceptable syringe | 50 μ L - 25mL | 同SMD02-1 |
| | Liquid volume | 1 μ L - 25mL | |
| | Precision | $\geq 0.5\%$ in the stroke of 30% - 100% | |
| | Accuracy | $\geq 99\%$ | |
| | Control resolution | 0.017% of syringe capacity | |
| | Min. volume increment | 0.017% of syringe capacity | |
| | Flow rates | 2.5 μ L/min to 1250ml/min | |
| | Communication | Bluetooth or RS232, Bluetooth wireless communication max. distance is 100 meters | |
| | Dimension (L×W×H) | 21.6×19.5×30.2 (cm) | |
| | Storage temperature | 20℃ - 65℃ | |
| | Operating temperature | 15℃ - 40℃ | |
| | Power supply | AC90 - 260V, 50/60Hz | |
| | Humidity | 10%-90%, (non-condensing) | |
| | Pipeline | PTFE tubing and borosilicate glass | |
| | Power consumption | <72W | <36W |
| | Control unit | 2 | 1 |
| | Max. Dilution proportion | 1:25000 | 1:50 |
| | Weight | 5.5kg | 3.5kg |