

Riffle Divider

Model: HGT- I / HGT- II / HGG- I / HGG- II

Model Difference:

Model	HGT-I	HGG-I	HGT-II	HGG-II
Material	Iron plate coating	stainless steel	Iron plate coating	stainless steel
Groove quantity	16		12	
Groove width	12.7mm		25.4mm	
Max sample classifying quantity	150 ~ 1200g		300 ~ 2000g	
Shape size	315 × 240 × 260mm		420 × 240 × 280mm	
Application range	middle and small seeds: rice, wheat, etc.		big seeds: bean, corn, etc.	
Specification	middle and small		big	



HGG series



HGT series