




NAVID KALA DARAN CO.

Name	Automatic Seed Counter
Model	SLY
Action	Microcomputer automatic control, press key are touches, complete automated operation. LED shows the hypothesis number and the actual number. Count speed can change, no noise, the precision is higher. The machine looks beautiful, shell is made of metal.
Technical specification	<ul style="list-style-type: none"> - Measuring range: 1-99999 Application scope: small sample: 0.7-4mm x 0.7-12mm big sample: 3-10mm x 3-12mm - Time for test: more than 1000/3minute - Precision: big sample: 2/1000 small sample: 4/1000 -Power: 220V 50HZ - Continuous working time: More than 5 hours
More	<p>SLY- A: For small sample such as rice, wheat, rape, vegetables, flowers and plants, tobacco.</p> <p>SLY- B: For big sample such as corn, bean, Sunflower seed</p> <p>SLY-C: Both big sample & small sample.</p>



Tel: 98-21-66486437/66486438 Fax: 98-21-66485558 Internet E-mail: info@navidkala.com

www.navidkala.com

